

Connecting Roger AudioHub to Roger inspiro

Reference Sheet



School: _____

Roger contact: _____


Date: _____

Teacher: _____

Phone: _____









Serial number: _____

Connecting

1. Plug Roger AudioHub into power source.
2. Hold Roger microphone within 4 inches (10 cm) of Roger AudioHub. Press **Connect**.
3. Connect Roger AudioHub to the audio source by plugging the audio cable into **Audio IN** on the back of the Roger AudioHub.
If plugging in only one source, use socket #1. Plug the other end of the cable into the audio source, such as MP3, TV, PC, etc.
4. Press and hold the Roger AudioHub's power button  for 2 seconds. Then start the audio source signal.
5. Volume should be adjusted first on the audio source's, and then on the Roger AudioHub.



Roger AudioHub indicator light

-  Solid red – Roger AudioHub is off
-  Solid blue – Standard operation mode
-  Blinking blue – Roger microphone is off or out of range
-  Solid green – Microphone is muted
-  Green to red – Volume is adjusting +/-
-  Blinking yellow – Software update
-  Blinking red – Call your Roger Contact
-  Alternating red/yellow – Cable is not correctly plugged in