# Phonak

# Guide pratique rapide

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Type de test	Outil	Référence
Questionnaires sur les attentes	Outil COAT Characteristic of Amplification Tool	Sandridge, S. et Newman, C. (2006). Improving the Efficiency and Accountability of the Hearing Aid Selection Process – Use of the COAT. <i>Audiology online</i> , consulté le 02/11/2019 sur <a href="https://www.audiologyonline.com/articles/improving-efficiency-and-accountability-hearing-995">https://www.audiologyonline.com/articles/improving-efficiency-and-accountability-hearing-995</a>
Questionnaires sur l'analyse des besoins	Échelle COSI Client-Orientated Scale of Improvement	Individualized, based on up to five user-nominated goals, categorized and with improvement subjectively rated. Dillon, H., James, A. et Ginis, J. (1997) Disponible sur <a href="https://www.nal.gov.au/products/downloadable-software/cosi-and-hauq/">https://www.nal.gov.au/products/downloadable-software/cosi-and-hauq/</a>
	Profil GHABP Glasgow Hearing Aid Benefit Profile (GHABP)	Gatehouse, S. (1999). Glasgow Hearing Aid Benefit Profile: Derivation and validation of client-centered outcome measures for hearing aid services. <i>Journal of the American Academy of Audiology</i> 10, 80-103.
	Échelle SSQ Speech, Spatial and Qualities of Hearing scale (SSQ-12)	Noble, W., Søgaard Jensen, N., Naylor, G., Bhullar, N. et Akeroyd, M. (2013). A short form of the Speech, Spatial and Qualities of Hearing scale suitable for clinical use: The SSQ12. <i>Int. J. Audiol.</i> 52(6), 409-412.
Mesures pour les accompagnants	Échelle SOS-HEAR Significant Other Scale for Hearing Disability	Scarinci, N., Worrall, L., Hickson, L. (2009) Test disponible sur <a href="https://shrs.uq.edu.au/communication-disability-centre">https://shrs.uq.edu.au/communication-disability-centre</a>
	Évaluation et solutions FOCAS Family-Oriented Communication Assessment and Solutions	Crowhen, D. et Turnbull, B. (2018) Test disponible sur <a href="https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/counseling_tools/documents/interactive_focas.pdf">https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/counseling_tools/documents/interactive_focas.pdf</a> Mode d'emploi disponible sur <a href="https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/counseling_tools/documents/how_to_use_focas.pdf">https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/counseling_tools/documents/how_to_use_focas.pdf</a>

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Outil	Référence
<b>Échelle COSI</b> Client-Orientated Scale of Improvement	Individualized, based on up to five user nominated goals, categorized and with improvement subjectively rated. Dillon, H., James, A. et Ginis, J. (1997) Disponible sur <a href="https://www.nal.gov.au/products/downloadable-software/cosi-and-hauq/">https://www.nal.gov.au/products/downloadable-software/cosi-and-hauq/</a>
<b>Représentation TELEGRAM</b> Telephone, Employment, Legislation, Entertainment, Groups, Recreation, Alarms, Members of the family	A graphical presentation of hearing needs that can be completed before and after any intervention and incorporates a broad range of situations. Thibodeau, L. (2004).  Test disponible sur <a href="https://www.utdallas.edu/hhlab/resources-and-publications/clinical-tools/">https://www.utdallas.edu/hhlab/resources-and-publications/clinical-tools/</a>
<b>Évaluation et solutions FOCAS</b> Family-Oriented Communication Assessment and Solutions	Crowhen, D. et Turnbull, B. (2018) Test disponible sur <a href="http://www.phonakpro.com">www.phonakpro.com</a> (voir liens ci-dessus, dans le « Tableau des outils utiles pour l'évaluation des besoins de communication »)

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## Références

- American Speech Language Hearing Association (2002). Guidelines for Fitting and Evaluation of FM Systems. *ASHA Desk Reference*. Disponible sur <https://www.asha.org/policy/GL2002-00010.htm>
- Bench, J., Kowal, A. et Bamford, J. (1979). The BKB (Bamford-Kowal-Bench) sentence lists for partially hearing children. *British Journal of Audiology* 13(3):108-112
- Boothroyd, A. (1968). Developments in Speech Audiometry. *British Journal of Audiology* 7(3):368-368
- Boothroyd, A., Hanin, L. et Hnath, T. (1985). A sentence test of speech perception: reliability, set equivalence, and short-term learning. *CUNY Academic Works*. Consulté en ligne le 02/11/2019, sur [https://academicworks.cuny.edu/cgi/viewcontent.cgi?article=1443&text=gc\\_pubs](https://academicworks.cuny.edu/cgi/viewcontent.cgi?article=1443&text=gc_pubs)
- British Society of Audiology (2018). Practice Guidance on the verification of hearing devices using probe microphone measurements. Disponible sur <http://www.thebsa.org.uk>
- British Association of Audiologists: (BAA Guideline) It is time to talk about Cochlear Implants. [https://www.baaudiology.org/files/8515/6267/2610/CI\\_BAA\\_Dickinson\\_FINAL\\_BAAtitle4.pdf](https://www.baaudiology.org/files/8515/6267/2610/CI_BAA_Dickinson_FINAL_BAAtitle4.pdf)
- Carlsson, P., Hjaldaahl J., Magnuson, A., Terneval, E., Eden, M., Skagerskarand, A. et Jonsson, R. (2014). Severe to profound hearing impairment: quality of life, psychosocial consequences and audiological rehabilitation. *Disability & Rehabilitation*, Early Online: 1-8
- Crowhen, D. et Turnbull, B. (2018). FOCAS: Family Oriented Communication Assessment and Solutions – A new holistic tool for performance hearing needs assessments. *Hearing Review*, Novembre : 20-26. <http://www.hearingreview.com/2018/10/focas-family-oriented-communication-assessment-solutions/>
- Dillon, H., James, A. et Ginis, J. (1997). Client Oriented Scale of Improvement (COSI) and its relationship to several other measures of benefit and satisfaction provided by hearing aids. *Journal of the American Academy of Audiology*. 8:27-43.
- Gatehouse, S. (1999). Glasgow Hearing Aid Benefit Profile: Derivation and validation of client-centered outcome measures for hearing aid services. *Journal of the American Academy of Audiology*. 10:80-103.
- Gopinath, B., Schneider, J., McMahon, C. M., Burlutsky, G., Leeder, S. R. et Mitchell, P. (2013). Dual sensory impairment in older adults increases the risk of mortality:

- a population-based study. *PloS One*, 8(3), e55054.  
DOI: 10.1371/journal.pone.0055054
- Grenness, C., Hickson, L., Laplante-Levesque, A., Davidson, B. (2014). Patient-centered care – A review for rehabilitative audiologists. *International Journal of Audiology*, 53(S), 60-67.
- Grimby, A. et Ringdah, A. (2000). Does Having a Job Improve the Quality of Life among Post-Lingually Deafened Swedish Adults with Severe-Profound Hearing Impairment?, *British Journal of Audiology*, 34:3, 187-195, DOI: 10.3109/03005364000000128
- Hallam, R.S., Jakes, S.C., Hinchcliffe, R. (1988). Cognitive variables in tinnitus annoyance. *Brit. J. Clin. Psychol.*; 27:213-222
- Hallam, R., Ashton, P., Sherbourne, K., Gailey, L. (2006). Acquired profound hearing loss: Mental health and other characteristics of a large sample, *International Journal of Audiology* 45(12), 715-723, DOI: 10.1080/14992020600957335
- Henry, J.A., Stewart, B.J., Abrams, H.B., Newman, C.W., Griest, S., Martin, W.H., Myers, P.J., Searchfield, G. (2014). Tinnitus Functional Index – Development and Clinical Application. *Audiology Today* 26(6), 40-48
- Newman, C.W., Jacobson, G.P., Spitzer, J.B. (1996). Development of the Tinnitus Handicap Inventory. *Arch.Otolaryngology* 122, 143-148
- NICE Guideline: Hearing loss in adults: assessment and management (2018). [www.nice.org.uk/guidance/ng98](http://www.nice.org.uk/guidance/ng98)
- Nilsson, M., Soli, S. et Sullivan, J. (1994). Development of the Hearing in Noise Test for the measurement of speech reception thresholds in quiet and in noise. *Journal of the Acoustical Society of America* 95(2), 1085-1099
- Niquette, P., Arcaroli, J., Revit, L., Parkinson, A., Staller, S., Skinner, M. et Killion, M. (2003). Development of the BKB-SIN Test. Article présenté lors de la réunion annuelle de l'*American Auditory Society* à Scottsdale, Arizona
- Noble, W., Søgaard Jensen, N., Naylor, G., Bhullar, N. et Akeroyd, M. (2013). A short form of the Speech, Spatial and Qualities of Hearing scale suitable for clinical use: The SSQ12. *International Journal of Audiology*. 52(6): 409-412.
- Peterson, G. et Lehiste, I. (1962). Revised CNC Lists for Auditory Tests. *Journal of Speech and Hearing Disorders*. 27(1)62
- Sandridge, S. et Newman, C. (2006). Improving the Efficiency and Accountability of the Hearing Aid Selection Process – Use of the COAT. *Audiology Online*, consulté le 02/11/2019 SUR <https://www.audiologyonline.com/articles/improving-efficiency-and-accountability-hearing-995>
- Scarinci, N., Worrall, L., Hickson, L. (2009). The Effect of Hearing Impairment in Older People on the Spouse: Development and Psychometric Testing of the Significant Other Scale for Hearing Disability (SOS-HEAR). *International Journal of Audiology* 48(10), 671-683.
- Schmitt, N., Winkler, A., Boretzki, M., Holube, I. (2016). A Phoneme Perception Test Method for High-Frequency Hearing Aid Fitting. *Journal of the American Academy of Audiology* 27(5), 367-379(13).  
<https://doi.org/10.3766/jaaa.15037>
- Spahr, A., Dorman, M., Litvak, L., Van Wie, S., Gifford, R., Loizou, P., Loiseau, L., Oakes, T. et Cook, S. (2012). Development and Validation of the AzBio Sentence Lists. *Ear and Hearing*. 33(1), 112-117
- Thibodeau, L. (2004). Maximizing communication via hearing assistance technology: Plotting beyond the audiogram! *Hearing Journal*, 57(11), 46-51.
- Turton, L. et Smith, P. (2013). Prevalence & characteristics of severe and profound hearing loss in adults in a UK National Health Service clinic. *International Journal of Audiology*, 52(2), 92-97.  
<https://doi.org/10.3109/14992027.2012.735376>
- Wagener, K.C., Brand T., Kollmeier B. (1999). Entwicklung und Evaluation eines Satztests für die deutsche Sprache Teil III: Evaluation des Oldenburger Satztests. *Z. Audiol.*, 38(3):86-95