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- Akeroyd M.A. & Summerfield A.Q. (2000). Integration of monaural and binaural evidence of vowel formants. *Journal of the Acoustical Society of America*, 107(6), 3394-3406.
- Balfour P.B. & Hawkins D.B. (1992). A comparison of sound quality judgments for monaural and binaural hearing aid processed stimuli. *Ear and Hearing*, 13(5), 331-339.
- Best V., Mason C.R. and Kidd G. (2011). Spatial release from masking in normally hearing and hearing-impaired listeners as a function of the temporal overlap of competing talkers. *Journal of the Acoustical Society of America* 129(3), 1616-1625.
- Boymans M., Goverts S.T., Kramer S.E., Festen J.M. & Dreschler W.A. (2008). A prospective multi-centre study of the benefits of bilateral hearing aids. *Ear and Hearing*, 29(6), 930-941.
- Bregman A.S. (1990). *Auditory Scene Analysis*. MIT Press: Cambridge, MA.
- Bronkhorst A.W. (2000). The cocktail party phenomenon: A review of research on speech intelligibility in multiple talker conditions. *Acustica*, 86, 117-128.
- Bronkhorst A.W. & Plomp R. (1988). The effect of head-induced interaural time and level differences on speech intelligibility in noise. *Journal of the Acoustical Society of America*, 83, 1508-1516.
- Brungart D. & Simpson B. (2002). The effects of spatial separation in distance on the informational and energetic masking of a nearby speech signal. *Journal of the Acoustical Society of America* 112, 664-676.
- Cherry E.C. (1953). Some experiments on the recognition of speech, with one and two ears. *Journal of the Acoustical Society of America*, 25, 975-979.
- Colburn H.S., Shinn-Cunningham B., Kidd G. Jr. & Durlach N. (2006). The perceptual consequences of binaural hearing. *International Journal of Audiology*, 45(Suppl 1), S34-44.
- Culling J.F., Hawley M.L. & Litovsky R.Y. (2004). The role of head-induced interaural time and level differences in the speech reception threshold for multiple interfering sound sources. *Journal of the Acoustical Society of America*, 116, 1057-1065.
- Dalton D.S., Cruickshanks K.J., Klein B.E.K., Klein R., Wiley T.L., et al. (2003). The impact of hearing loss on quality of life in older adults. *The Gerontologist*, 43(5), 661-668.
- Day G., Browning G. & Gatehouse S. (1988). Benefit from binaural hearing aids in individuals with a severe hearing impairment. *British Journal of Audiology*, 23, 273-277.
- Dillon A. (2001). Beyond usability: process, outcome and affect in human-computer interactions. *Canadian Journal of Library and Information Science*, 26(4), 57-69.
- Duquesnoy A.J. (1983). The intelligibility of sentences in quiet and in noise in aged listeners. *Journal of the Acoustical Society of America*, 74, 1136-1144.
- Durlach N.I., Mason C.R., Kidd G. Jr., Arbogast T.L., Colburn H.S. & Shinn-Cunningham B.G. (2003). Note on informational masking. *Journal of the Acoustical Society of America*, 113, 2984-2987.
- French N.R. & Steinberg J.C. (1947). Factors governing the intelligibility of speech sounds. *Journal of the Acoustical Society of America*, 19, 90-119.
- Haggard M. & Hall J. (1982). Forms of binaural summation and the implications of individual variability for binaural hearing aids. *Scandinavian Audiology Supplementum* (15), 47-63.
- Hawkins D.B., Walden B.E., Montgomery A. & Prosek R.A. (1987). Description and validation of an LDL procedure designed to select SSPL90. *Ear and Hearing*, 8(3), 162-169.
- Hornsby B.W., Ricketts T. A. & Johnson E. E. (2006). The effects of speech and speechlike maskers on unaided and aided speech recognition in persons with hearing loss. *Journal of the American Academy of Audiology*, 17, 432-447.
- Keys J.W. (1947). Binaural versus monaural hearing. *Journal of the Acoustical Society of America*, 19(4), 629-631.
- Kobler S. & Rosenhall U. (2002). Horizontal localization and speech intelligibility with bilateral and unilateral hearing aid amplification. *International Journal of Audiology*, 41, 395-400.
- Kochkin S. & Kuk F. (1997). The binaural advantage: evidence from subjective benefits and customer satisfaction data. *The Hearing Review*, 4(4), 29-32.
- Kochkin S. (2000). MarkeTrak V: Consumer Satisfaction Revisited. *The Hearing Journal*, 53(1), 38-55.
- Kreikemeier S., Margolf-Hackl S., Raether J., Fichtl E. & Kiessling J. (2012). Vergleichende Evaluation unterschiedlicher Hörgeräte-Richtmikrofontechnologien bei hochgradig Schwerhörigen. *Zeitschrift für Audiologie, Supplement zur 15. Jahrestagung der deutschen Gesellschaft für Audiologie*.
- McArdle R., Killion M., Mennite M. & Chisolm T. (2012). Are two ears not better than one? *Journal of the American Academy of Audiology*, 23, 171-181.
- Nabelek A.K. & Pickett J.M. (1974). Reception of consonants in a classroom as affected by monaural and binaural listening, noise, reverberation, and hearing aids. *Journal of the Acoustical Society of America*, 56, 628-639.
- Noble W. & Gatehouse S. (2006). Effects of bilateral versus unilateral hearing aid fitting on abilities measured by the Speech, Spatial, and Qualities of Hearing Scale (SSQ). *International Journal of Audiology*, 45, 172-181.
- Nyffeler M. (2010). DuoPhone. *Field Study News*, February.
- O'Neil J.N., Connelly C.J., Limb C.J. & Ryugo D.K. (2011). Synaptic morphology and the influence of auditory experience. *Hearing Research*, 279(102), 118-130.
- Oxenham A.J. & Plack C.J. (1997). A behavioral measure of basilar membrane nonlinearity in listeners with normal and impaired hearing. *Journal of the Acoustical Society of America*, 101, 3666-3675.
- Peissig J. & Kollmeier B. (1997). Directivity of binaural noise reduction in spatial multiple noise-source arrangements for normal and impaired listeners, *Journal of the Acoustical Society of America*, 101, 1660-1670.
- Pichora-Fuller M.K. & Singh G. (2006). Effects of age on auditory and cognitive processing: implications for hearing aid fitting and audiological rehabilitation. *Trends in Amplification*, 10, 29-59.
- Picou E.M. & Ricketts T.A. (2011). Comparison of wireless and acoustic hearing aid-based telephone listening strategies. *Ear and Hearing*, 32(2), 209-220.
- Reynolds G.S. & Stevens S.S. (1960). Binaural summation of loudness. *The Journal of the Acoustical Society of America*, 32(10), 1337-1344.
- Singh G., Pichora-Fuller M.K. & Schneider B.A. (2008). The effect of age on auditory spatial attention in conditions of real and simulated spatial separation. *Journal of the Acoustical Society of America*, 124, 1294-1305.
- Sivonen V.P. (2011). Binaural directivity patterns for normal and aided human hearing. *Ear & Hearing*, 32, 674-677.
- Shaw E.A.G. (1974). Transformation of sound pressure level from the free field to the eardrum in the horizontal plane. *Journal of the Acoustical Society of America*, 56, 1848-1860.
- Yueh B., Shapiro N., MacLean C.H., et al. (2003, April 16). Screening and management of adult hearing loss in primary care. *Journal of the American Medical Association*, 289 (15), 1976-1985.
- Zurek P. M. (1980). The precedence effect and its possible role in the avoidance of interaural ambiguities. *Journal of the Acoustical Society of America*, 67, 953-964.