Phonak Field Study News.

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- Abrams, H. B. & Kihm, J. (2015). An introduction to MarkeTrak IX: A New Baseline for the Hearing Aid Market. *Hearing Review*, 22(6).
- Appleton, J. (2022). What Is Important to Your Hearing Aid Clients...and Are They Satisfied? *Hearing Review*. 29 (6).
- Appleton, J. (2020) AutoSense OS 4.0 – significantly less listening effort and preferred for speech intelligibility. Phonak Field Study News retrieved from www.phonakpro.com/evidence, accessed October 31st, 2022.
- Appleton, J. & König, G. (2014). Improvements in speech intelligibility and subjective benefit with binaural beamformer technology. *Hearing Review*, 21(11), 40–42.
- Derleth, P., Georganti, E., Latzel, M., Courtois, G., Hofbauer, M., Raether, J., & Kuehnelt, V. (2021). Binaural Signal Processing in Hearing Aids. *Seminars in Hearing*, 42, 206 – 223.
- Hagerman, B. & Olofsson, A. (2004). A method to measure the effect of noise reduction algorithms using simultaneous speech and noise. *Acta Acustica united with Acustica*. 90, 356–361.
- Hayes, D. (2019). Speech detection by direction. Unitron White Paper. Retrieved from: https://www.unitron.com/au/en_au/learn/speech-detection-by-direction0.html, accessed October 31st, 2022.
- Hornsby, B.W.Y. (2013). The effects of hearing aid use on listening effort and mental fatigue associated with sustained speech processing demands. *Ear and Hearing*, 34(5), 523–534.
- Jensen N. S., Neher T., Laugesen S., Johannesson, R. B. & Kragelind, L. (2013). Laboratory and field study of the potential benefits of pinna cue-preserving hearing aids. *Trends in Amplification*. 17 (3/4). 171–188.
- Jespersen, C. T., Kirkwood, B. C., Groth, J. (2021). Increasing the effectiveness of hearing aid directional microphones. *Seminars in Hearing*. 42: 224–236.
- Keidser, G., O'Brian, A., Hain, J., McLelland, M., & Yeend, I. (2009). The effect of frequency-dependent microphone directionality on horizontal localization performance in hearing-aid users. *International Journal of Audiology*, 48(11), 789–803.
- Killion, M. C. (1997). The SIN report: Circuits haven't solved the hearing-in-noise problem. *Hearing Journal*, 50(10), 28–32.
- Latzel, M., & Appleton-Huber, J. (2015). StereoZoom – Adaptive behaviour improves speech intelligibility, sound quality and suppression of noise. *Field Study News*. Retrieved from www.phonakpro.com/evidence, accessed October 31st, 2022.
- Luts, H., Eneman, K., Wouters, J., Schulte, M., Vormann, M., Buechler, M., Dillier, N., Houben, R., Dreschler, W. A., Froehlich, M., Puder, H., Grimm, G., Hohmann, V., Leijon, A., Lombard, A., Mauler, D., & Spriet, A. (2010). Multicenter evaluation of signal enhancement algorithms for hearing aids. *The Journal of the Acoustical Society of America*, 127(3), 1491–1505.
- Krüger, M., Schulte, M., Brand, T., & Holube, I. (2017). Development of an adaptive scaling method for subjective listening effort. *The Journal of the Acoustical Society of America*, 141(6), 4680–4693.
- Pichora-Fuller, M.K., & Singh, G. (2006). Effects of Age on Auditory and Cognitive Processing: Implications for Hearing Aid Fitting and Audiologic Rehabilitation. *Trends in Amplification*, 10(1), 29–59.
- Picou, E.M., Ricketts, T.A., & Hornsby, B.W.Y. (2013). How hearing aids, background noise, and visual cues influence objective listening effort. *Ear and Hearing*, 34(5), e52–e64.
- Picou, E., Aspell, E., & Ricketts, T. (2014). Potential benefits and limitations of three types of directional processing in hearing aids. *Ear and Hearing*, 35(3), 339–52.
- Raether, J. (2005). Real Ear Sound – A simulation of the pinna effect optimizes sound localization also with open fittings. *Phonak Field Study News*. Retrieved from www.phonakpro.com/evidence, accessed October 31, 2022.
- Schulte, M., Meis, M., Krüger, M., Latzel, M., Appleton-Huber, J. (2018). Significant increase in the amount of social interaction when using StereoZoom. *Phonak Field Study News*. Retrieved from: www.phonakpro.com/evidence, accessed October 31st, 2022.
- Stewart, E. Rakita, L., & Drexler, J. (2019). Phonak Compendium: StereoZoom Part 1: The benefit of wirelessly connected narrow directionality in Phonak hearing aids for speech intelligibility. Retrieved from: www.phonakpro.com/evidence, accessed October 31st, 2022.
- Xu, J., & Han, W. (2014). Improvement of Adult BTE Hearing Aid Wearers' Front/Back Localization Performance Using Digital Pinna-Cue Preserving Technologies: An Evidence-Based Review. *Korean Journal of Audiology*, 18, 97–104.

Walden B. E., Surr R. K., Cord M.T., & Dyrland O. (2004). Predicting hearing aid microphone preference in everyday listening. *J Am Acad Audiol.* 15(5): 365-96.

Wagener, K., Kühnel, V., and Kollmeier, B. (1999). Development and evaluation of a German sentence test I: Design of the Oldenburg sentence test. *Z. Audiol.* 38, 4-15.

Wagener, K. C., & Brand, T. (2005). Sentence intelligibility in noise for listeners with normal hearing and hearing impairment: influence of measurement procedure and masking parameters. *International journal of audiology*, 44(3), 144-156.

Winneke, A., Schulte, M., Vormann, M. & Latzel, M. (2020). Effect of directional microphone technology in hearing aids on neural correlates of listening and memory effort: an electroencephalographic study. *Trends in Hearing.* 24:1-16. Retrieved from www.phonakpro.com/evidence, accessed October 31st, 2022.

Woodward, J., Kühnel, V., & Latzel, M. (2022). Innovative speech understanding solutions: Lumity with SmartSpeech™ Technology. Phonak Insight. Retrieved from: www.phonapro.com/evidence, accessed October 15th, 2022.

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028-2475-02/V1.00/2019-08/kt © 2022 Sonova AG All rights reserved