# Field Study News

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---

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Chung, K., 2004. Challenges and Recent Developments in Hearing Aids. Part I. Speech Understanding in noise, microphone technologies and noise reduction algorithms. *Trends in amplification*, 8(3), p. 83-124.

Hamacher, V., Eggers, J., Fischer, E., Kornagel, U., Puder, H., Rass, U., 2005. Signal Processing in High-End Hearing Aids: State of the Art, Challenges, and Future Trends EURASIP. *Journal of Applied Signal Processing*, 18 p. 2915-2929

Latzel, M., 2013. Concepts for Binaural Processing in Hearing Aids. *Hearing Review*, 20(4), p. 34

Nyffeler, M., 2010. StereoZoom - Improvements with directional microphones. Field Study News. Phonak AG: 2010

Picou, E. M., Aspell, E., Ricketts, T. A., 2014. Potential benefits and limitations of directional processing in hearing aids. *Ear and Hearing*, 35(3), p. 339-352

Ricketts, T. A., Henry, P., 2002. Evaluation of an adaptive, directional-microphone hearing aid. *International Journal of Audiology*, 41 p. 100-112

Ricketts, T. A., 2006. Directional hearing aid benefit in listeners with severe hearing loss. *International Journal of Audiology*, 45, p. 190-197

Stuermann, B., 2011. StereoZoom - Improved speech understanding even with open fittings. Field Study News. Phonak AG: 2011

Wouters, J., Litierère, L., van Wieringen, A., 1999. Speech intelligibility in noisy environments with one- and two-microphone hearing aids. *Audiology*, 38 p. 91-98

Wouters, J., Vanden Berghe, J., Maj, J.-B., 2002. Adaptive noise suppression for a dual-microphone hearing aid. *International Journal of Audiology*, 41 p. 401-407