Phonak

Field Study News

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Abrams, H. B., & Kihm, J. (2015). An Introduction to MarkeTrak IX: A New Baseline for the Hearing Aid Market. *Hearing Review*, 22(6).

Helbling, T., Vormann, M., & Wagener, K. (2013). Speech priority noise reduction: A digital noise reduction preference study. *Hearing Review*, 10(11), 34–43.

Hochmuth, S., Jürgens, T., Brand, T., & Kollmeier, B. (2015). Talker- and language-specific effects on speech intelligibility in noise assessed with bilingual talkers: Which language is more robust against noise and reverberation? *International Journal of Audiology*, 54(2), 23–34.

Kalbe E1, Kessler J, Calabrese P, Smith R, Passmore AP, Brand M, Bullock R. (2004). DemTect: a new, sensitive cognitive screening test to support the diagnosis of mild cognitive impairment and early dementia. *International Journal of Geriatric Psychiatry*, 19(2), 136–143.

Killion, M. C. (1997). The SIN report: Circuits haven't solved the hearing-in-noise problem. *Hearing Journal*, 50(10), 28–32.

Meis, M., Krueger, M., Gablenz, P. v., Holube, I., Gebhard, M., Latzel, M., & Paluch, R. (2018). Development and application of an annotation procedure to assess the impact of hearing aid amplification on interpersonal communication behavior. *Trends in Hearing*, 22, 1–17.

Paluch, R., Krüger, M., Hendrikse, M., Grimm, G., Hohmann, V., & Meis, M. (2017). Ethnographic research: The interrelation of spatial awareness, everyday life/laboratory environments and effects of hearing aids. *Proceedings of the International Symposium on Auditory and Audiological Research (Proc. ISAAR)*, Vol. 6: Adaptive Processes in Hearing, August 2017, Nyborg, Denmark. Edited by Santurette, S., Dau, T., Dalsgaard, J. C., Tranebjærg, L., Andersen, T., & Poulsen, T. The Danavox Jubilee Foundation, 2017. ©

Shinn-Cunningham, B. (2008). Objective-based auditory and visual attention. *Trends in Cognitive Sciences*, 12(5), 182–186.

Wagener, K. C., Hochmuth, S., Ahrlich, M., Zokoll, M. A., & Kollmeier, B. (2014). Der weibliche Oldenburger Satztest 17th annual conference of the DGA Oldenburg, CD-Rom.

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Authors and investigators

External principle investigators

Since 2004, Michael Schulte has been with Hörzentrum Oldenburg GmbH, Germany, where he has been responsible for audiological studies in publicly funded projects as well as in cooperation with the industry. In 2002, he received his Ph.D. from the Biomagnetism

Centre at the Institute of Experimental Audiology, University of Münster, Germany. From 2002 to 2003, he worked as a postdoc at the F.C. Donders Centre for Cognitive Neuroimaging, Nijmegen, Netherlands. Michael Schulte's research interest is in the evaluation of hearing systems with a special focus on listening effort.

Dr. Matthias Vormann received his diploma in physics in 1995 and his Ph.D. in 2011 from University of Oldenburg. Since 2005 he has been with Hörzentrum Oldenburg GmbH, Germany, where he works mainly in industrial research projects for

hearing aid manufacturers. His interest is about new methods for measuring the subjective and objective benefit of hearing aids.

Jan Heeren studied Physics at the University of Oldenburg, Germany, and graduated in the Medical Physics group in 2014. From 2012, he worked on several projects in the field of hearing aid evaluation and virtual acoustics at the university and the Hörzentrum Oldenburg. In

2016, he started in the R&D department at HörTech GmbH, Oldenburg, working on hearing aid evaluation methods. Apart from his scientific activities, he has conducted more than 500 events as a free-lancing audio engineer since 2008.

Study coordinator

Matthias Latzel studied electrical engineering in Bochum and Vienna in 1995. After completing his Ph.D. in 2001, he carried out his Postdoc from 2002 to 2004 in the Department of Audiology at Giessen University. He was the head of the Audiology department at Phonak

Germany from 2011. Since 2012 he has been working as the Clinical Research Manager for Phonak AG, Switzerland

Author

Jennifer Appleton-Huber received her M.Sc. in Audiology from the University of Manchester in 2004. Until 2013, she worked as an Audiological Scientist mainly in the UK and Switzerland, where she worked with adults and pediatrics,

in the areas of hearing aids and cochlear implants. Her current role is Technical Editorial Manager at Phonak Headquarters.