Pediatric Focus 3

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HEARING AID FITTING AND VERIFICATION CHECKLIST

CHILD: _____

DATE: _____

NOTES: _____

DESCRIPTION OF EAR CANAL ACOUSTICS

Transducer used to assess hearing thresholds: insert earphones + personal earmold
 insert earphones + foam-tip
 Other: _____

RECD for verification: new previously measured
RECD Coupler: HA-1 HA-2 0.4cc WRECD
RECD Coupling type: foam-tip personal earmold

If predicted RECD used, provide reason: _____

ELECTROACOUSTIC VERIFICATION OF FIT-TO-TARGETS AND SII VALUES

Soft level speech (55 dB SPL) R ear within ± 5 dB of DSL targets over targets under targets
 SII within normative range
L ear within ± 5 dB of DSL targets over targets under targets
 SII within normative range

Average level speech (65 dB SPL) R ear within ± 5 dB of DSL targets over targets under targets
 SII within normative range
L ear within ± 5 dB of DSL targets over targets under targets
 SII within normative range

Maximum power output (MPO) R ear within ± 5 dB of DSL targets over targets under targets
 SII within normative range
L ear within ± 5 dB of DSL targets over targets under targets
 SII within normative range

CONSIDERATION OF ADVANCED AMPLIFICATION TECHNOLOGIES

Noise management Candidate?: yes no Feature Enabled?: yes no
Verification documented? yes no _____ dB of noise reduction

Frequency lowering Candidate?: yes no Feature Enabled?: yes no
Verification documented? yes no

Remote microphone Candidate?: yes no Feature Enabled?: yes no
Verification documented? yes no

Feedback suppression Candidate?: yes no Feature Enabled?: yes no
 considered status of earmold(s) New earmold(s) required? yes no

Directional microphone Microphone mode selected: pinna matched fixed adaptive

Data-logging Feature Enabled?: yes no hrs/day of use: _____

Excerpt from Ontario Infant Hearing Program Amplification Protocol | 2019.01

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American Academy of Audiology (2011). American Academy of Audiology Clinical Practice Guidelines: Remote Microphone Hearing Assistance Technologies for Children and Youth from Birth to 21 Years. Retrieved from https://audiology-web.s3.amazonaws.com/migrated/HAT_Guidelines_Supplement_A.pdf_53996ef7758497.54419000.pdf

American Academy of Audiology (2011). American Academy of Audiology Clinical Practice Guidelines: Remote Microphone Hearing Assistance Technologies for Children and Youth from Birth to 21 Years. Retrieved from

https://audiology-web.s3.amazonaws.com/migrated/HAT_Guidelines_Supplement_A.pdf_53996ef7758497.54419000.pdf

American Academy of Audiology (2013). American Academy of Audiology Clinical Practice Guidelines on Pediatric Amplification. Retrieved from <http://www.audiology.org/resources/documentlibrary/Documents/PediatricAmplificationGuidelines.pdf>

American Academy of Audiology. (2020). *Clinical Guidance Document: Assessment of Hearing in Infants and Young Children* (p. 56) [Clinical Guidance Document]. American Academy of Audiology. https://www.audiology.org/sites/default/files/publications/resources/Clin%20Guid%20Doc_Assess_Hear_Infants_Children_1.23.20.pdf

Bagatto, M. (2020). Audiological Considerations for Managing Mild Bilateral or Unilateral Hearing Loss in Infants and Young Children. *Language, Speech, and Hearing Services in Schools, 51*(1), 68–73. https://doi.org/10.1044/2019_LSHSS-OCHL-19-0025

Bagatto, M., Moodie, S., Brown, C., Malandrino, A., Richert, F., Clench, D., & Scollie, S. (2016). Prescribing and verifying hearing aids applying the American Academy of Audiology Pediatric Amplification guideline: Protocols and outcomes from the Ontario Infant hearing program. *Journal of the American Academy of Audiology, 27*(3). <https://doi.org/10.3766/jaaa.15051>

Bagatto, M., Moodie, S., Scollie, S., Seewald, R., Moodie, S., Pumford, J., & Liu, K. P. R. (2005). Clinical protocols for hearing instrument fitting in the desired sensation level method. *Trends in Amplification, 9*(4). <https://doi.org/10.1177/108471380500900404>

Bagatto, M. P., Moodie, S. T., Malandrino, A. C., Richert, F. M., Clench, D. A., & Scollie, S. D. (2011). The University of Western Ontario pediatric audiological monitoring protocol (UWO PedAMP). *Trends in Amplification, 15*(1). <https://doi.org/10.1177/1084713811420304>

Benítez-Barrera CR, Angley GP, Tharpe AM. (2017). “Remote Microphone Use at Home: Impact on Caregiver Talk”. *Journal of Speech, Hearing, Language Research, 61*(2): 399-409. https://doi.org/10.1044/2017_JSLHR-H-17-0168.

Benítez-Barrera, C.R., Thompson, E.C., Angley, G.P., & Tharpe, A.M. (2019). Remote Microphone System Use at Home: Impact on Child-Directed Speech. *Journal of Speech, Hearing, Language Research, 62*: 2002-2008. https://doi.org/10.1044/2019_JSLHR-H-18-0325

British Society of Audiology. (2014). *Recommended Procedure: Visual Reinforcement Audiometry* (Recommended Procedure No. OD104-37; pp. 1–27). <https://www.thebsa.org.uk/wp-content/uploads/2014/06/OD104-37-Recommended-Procedure-Visual-Reinforcement-Audiometry-2014-1.pdf>

British Society of Audiology. (2019). *Auditory Brainstem Response (ABR) Testing in Babies* (Recommended Procedure OD-104-81; pp. 1–58). British Society of Audiology. <https://www.thebsa.org.uk/wp-content/uploads/2019/06/OD104-81-Recommended-Procedure-for-ABR-Testing-in-Babies.pdf>

Ching, T.Y.C., Quar, T.K., Johnson, E.E., Newall, P., Sharma, M. (2015). Comparing NAL-NL1 and DSL v5 in Hearing Aids Fit to Children with Severe or Profound Hearing Loss: Goodness of Fit-to-Targets, Impacts on Predicted Loudness and Speech Intelligibility. *Journal of the American Academy of Audiology, 26*(3): 260-274.

Ching, T.Y.C., Zhang, V.W., Johnson, E.E., Van Buynder, P., Hou, S., Burns, L., Button, L., Fynn, C., McGhie, K. (2018). Hearing aid fitting and developmental outcomes of children fit according to either the NAL or DSL prescription: fit-to-target, audibility, speech and language abilities. *International Journal of Audiology, 27*(sup2): S41-S54.

Cruckley, J., & Scollie, S. D. (2012). Children’s speech recognition and loudness perception with the Desired Sensation Level v5 Quiet and Noise prescriptions. *American Journal of Audiology, 21*(2). [https://doi.org/10.1044/1059-0889\(2012/12-0002\)](https://doi.org/10.1044/1059-0889(2012/12-0002))

- Crukley, J., Scollie, S., & Parsa, V. (2011). An Exploration of Non-Quiet Listening at School. *Journal of Educational Audiology*, 17, 23–35.
- Curran, M., Walker, E. A., Roush, P., & Spratford, M. (2019). Using propensity score matching to address clinical questions: the impact of remote microphone systems on language outcomes in children who are hard of hearing. *Journal of Speech, Language, and Hearing Research*, 62(3), 564–576. https://doi.org/10.1044/2018_JSLHR-L-ASTM-18-0238
- Davidson Lisa S., & Skinner Margaret W. (2006). Audibility and Speech Perception of Children Using Wide Dynamic Range Compression Hearing Aids. *American Journal of Audiology*, 15(2), 141–153. [https://doi.org/10.1044/1059-0889\(2006/018\)](https://doi.org/10.1044/1059-0889(2006/018))
- EUHA. (2015). *Zusammenfassung EUHA-MPO-Signal und Einstellung der MPO* (pp. 1–2). <http://www.euha.org/assets/Uploads/Leitlinien/Expertenkreis-04-Hoerakustik/EUHA-Leitlinie-04-01-Ergaenzung-1.pdf>
- Feilner, M., Rich, S., & Jones, C. (2016). *Automatic and directional for kids* (pp. 1–5) [Phonak Insight]. https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/evidence/white_paper/documents/technical_paper/insight_automatic_and_directional_for_kids_028-1499.pdf
- Folkeard, P., Bagatto, M., & Scollie, S. (2020). Evaluation of Hearing Aid Manufacturers' Software-Derived Fittings to DSL v5.0 Pediatric Targets. *Journal of the American Academy of Audiology*, 31(5), 354–362. <https://doi.org/10.3766/jaaa.19057>
- Folkeard, P., Pumford, J., Abbasalipour, P., Willis, N., & Scollie, S. (2018). A Comparison of Automated Real-Ear and Traditional Hearing Aid Fitting Methods—Audioscan. *Hearing Review*, 25(11), 28–32.
- Glista, D., Hawkins, M., Bohnert, A., Rehmann, J., Wolfe, J., & Scollie, S. (2017). The effect of adaptive nonlinear frequency compression on phoneme perception. *American Journal of Audiology*, 26(4). https://doi.org/10.1044/2017_AJA-17-0023
- Glista, D., Hawkins, M., Scollie, S., Wolfe, J., Bohnert, A., & Rehmann, J. (2016). *Pediatric verification for SoundRecover2* (Best Practice Protocol No. 028-1528-03/V1.00/2016-05; pp. 1–8). Phonak.
- Glista, D., Hawkins, M., Scollie, S., Wolfe, J., Bohnert, A., & Rehmann, J. (2016). Best practice protocol: Pediatric verification for SoundRecover2. Phonak, Accessed from: https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/products_solutions/pediatrics/documents/best_practice_protocols/best_practice_protocol_sound_recover2_pediatric_verification.pdf.
- Glista, D., Scollie, S., Moodie, S., Easwar, V., and Network of Pediatric Audiologists (2014). The Ling 6 (HL) test: Typical pediatric performance data and clinical use evaluation, *Journal of the American Academy of Audiology*, 25(10): 1008-1021
- Glista, D., Scollie, S., & Sulkers, J. (2012). Perceptual acclimatization post nonlinear frequency compression hearing aid fitting in older children. *Journal of Speech, Language, and Hearing Research*, 55(6). <https://doi.org/10.1044/1092-4388>
- Gustafson, S. J., & Pittman, A. L. (2011). Sentence perception in listening conditions having similar speech intelligibility indices. *International Journal of Audiology*, 50(1), 34–40. <https://doi.org/10.3109/14992027.2010.521198>
- Gustafson, S. J., Ricketts, T. A., & Tharpe, A. M. (2017). Hearing Technology Use and Management in School-Age Children: Reports from Data Logs, Parents, and Teachers. *Journal of the American Academy of Audiology*, 28(10), 883–892. <https://doi.org/10.3766/jaaa.16042>
- Hands and Voices Inc. (2019). Guide for parent-professional partnerships in the hearing healthcare setting. Phonak.

https://www.phonak.com/content/dam/phonakpro/gc_hq/en/resources/counseling_tools/documents/Brochure_BtBtC_Guide_parents_and_professionals_Hand_and_Voices_210x280_EN_V1.00_028-1956-02.pdf

Hodgetts, W. E., & Scollie, S. D. (2017). DSL prescriptive targets for bone conduction devices: adaptation and comparison to clinical fittings. *International Journal of Audiology*, 56(7).
<https://doi.org/10.1080/14992027.2017.1302605>

Joint Committee on Infant Hearing. (2019). Year 2019 Position Statement: Principles and Guidelines for Early Detection and Intervention Programs. *The Journal of Early Hearing Detection and Intervention*, 4(2): 1-44.
<https://digitalcommons.usu.edu/cgi/viewcontent.cgi?article=1104&context=jehdi>

Keidser, G., Dillon, H. R., Flax, M., Ching, T., & Brewer, S. (2011). The NAL-NL2 prescription procedure. *Audiology Research*, 1(1S), 1–3. <https://doi.org/10.4081/audiores.2011.e24>

Latzel, M. Denys, S., Anderson, S., Francart, T., Wouters, J. & Appleton-Huber, J. (2017). An integrated REM system with proven accuracy and reliability. *The Hearing Review*, October.
<http://www.hearingreview.com/2017/09/integrated-rem-system-proven-accuracy-reliability/>

Lewis, D., & Bagatto, M. (2017). Considering directional microphone use in pediatric hearing aid fittings. *Pediatric Focus 1*, Phonak. Accessed from:
https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/evidence/white_paper/documents/technical_paper/focus_btb_pediatric_directional_microphone_use_in_pediatric_ha_fitting.pdf.

Marriage, J. E., Moore, B. C. J., Stone, M. A., & Baer, T. (2005). Effects of Three Amplification Strategies on Speech Perception by Children With Severe and Profound Hearing Loss. *Ear and Hearing*, 26(1), 35–47.

McCreery, R. W., Bentler, R. A., & Roush, P. A. (2013). Characteristics of hearing aid fittings in infants and young children. *Ear and Hearing*, 34(6), 701–710. <https://doi.org/10.1097/AUD.0b013e31828f1033>

McCreery, R., Brennan, M., Walker, E.A., & Spratford, M. (2017). Perceptual implications of level- and frequency-specific deviations from hearing aid prescription in children. *Journal of the American Academy of Audiology*, 28(9):861-875.

McCreery, R., Walker, E., Spratford, M., Bentler, R., Holte, L., Roush, P., ... Moeller, M. P. (2015). Longitudinal predictors of aided speech audibility in infants and children. *Ear and Hearing*, 36, 24S–37S.
<https://doi.org/10.1097/AUD.0000000000000211>

McCreery, R. W., Walker, E. A., Stiles, D. J., Spratford, M., Oleson, J. J., & Lewis, D. E. (2020). Audibility-Based Hearing Aid Fitting Criteria for Children With Mild Bilateral Hearing Loss. *Language, Speech, and Hearing Services in Schools*, 51(1), 55–67. https://doi.org/10.1044/2019_LSHSS-OCHL-19-0021

Moeller, M. P., & Tomblin, J. B. (2015). Epilogue: Conclusions and Implications for Research and Practice. *Ear and Hearing*, 36, 92S-98S. <https://doi.org/10.1097/AUD.0000000000000214>

Moodie, S., Pietrobon, J., Rall, E., Lindley, G., Eiten, L., Gordey, D., ... Scollie, S. (2016a). Using the real-ear-to-coupler difference within the American Academy of audiology pediatric amplification guideline: Protocols for applying and predicting earmold RECDs. *Journal of the American Academy of Audiology*, 27(3).
<https://doi.org/10.3766/jaaa.15086>

Moodie, S., Rall, E., Eiten, L., Lindley, G., Gordey, D., Davidson, L., ... Scollie, S. (2016b). Pediatric audiology in North America: Current clinical practice and how it relates to the American Academy of Audiology pediatric amplification guideline. *Journal of the American Academy of Audiology*, 27(3).
<https://doi.org/10.3766/jaaa.15064>

Moodie, S. T. F., Scollie, S. D., Bagatto, M. P., & Keene, K. (2017). Fit-to-targets for the desired sensation level version 5.0a hearing aid prescription method for children. *American Journal of Audiology*, 26(3).
https://doi.org/10.1044/2017_AJA-16-0054

Moodie, S. T., & Sindrey, D. (n.d.). *Hear On Videos*— Western University. Hear On Videos. Retrieved July 17, 2020, from <https://www.uwo.ca/nca/fcei/hearon/index.html>

Ontario Infant Hearing Program. (2018). *Protocol for Auditory Brainstem Response – Based Audiological Assessment (ABRA)* (M. Bagatto, V. Easwar, R. El-Naji, M. Hyde, V. Martin, M. Pigeon, D. Purcell, S. Scollie, & J. Witte, Eds.; pp. 1–69).

https://www.uwo.ca/nca/pdfs/clinical_protocols/2018.01%20ABRA%20Protocol_Oct%2031.pdf

Ontario Infant Hearing Program. (2019a). *Audiometric Assessment for Children Aged 6 to 60 months* (S. Scollie, M. Pigeon, M. Bagatto, J. Witte, & A. Malandrino, Eds.; pp. 1–44). Ontario Ministry of Children, Community, and Social Services, Infant Hearing Program.

https://www.uwo.ca/nca/pdfs/clinical_protocols/IHP_CBA%20Protocol_2019.01.pdf

Ontario Infant Hearing Program. (2019b). *Protocol for the Provision of Amplification* (M. Bagatto & S. Scollie, Eds.; 2019.01, pp. 1–98). Ontario Ministry of Children, Community, and Social Services, Infant Hearing Program.

https://www.uwo.ca/nca/pdfs/clinical_protocols/IHP_Amplification%20Protocol_2019.01.pda

Phonak. (n.d.-a). *Child Hearing Assessments*. PhonakPro. Retrieved July 17, 2020, from

<https://www.phonakpro.com/ca/en/resources/counseling-tools/pediatric/child-hearing-assessments.html>

Phonak. (n.d.-b). *Hearing Aids for Children: Frequently asked questions*. Phonak. Retrieved July 17, 2020, from

<https://www.phonak.com/ca/en/support/children-and-parents/children-faq/hearing-aid-faq.html>

Phonak. (2016). *TargetMatch Fitting Guide* (V1.00/2016-03; p. 4). Phonak AG.

https://www.phonakpro.com/content/dam/phonakpro/gc_us/en/products_solutions/hearing_aid/documents/venture/TargetMatch_Fitting_Guide_Target.pdf

Phonak AG. (2013). *Phonak Insight. Junior mode—The latest developments in Phonak Target™ Junior mode* (Phonak Insight 6, pp. 1–6).

https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/evidence/white_paper/documents/technical_paper/Insight_Junior_Mode_028-0983.pdf

Pittman, A. L., Pederson, A. J., & Rash, M. A. (2014). Effects of Fast, Slow, and Adaptive Amplitude Compression on Children's and Adults' Perception of Meaningful Acoustic Information. *Journal of the American Academy of Audiology*, 25(9), 834–847. <https://doi.org/10.3766/jaaa.25.9.6>

Pittman, A. L., & Stelmachowicz, P. G. (2003). Hearing loss in children and adults: audiometric configuration, asymmetry, and progression. *Ear and Hearing*, 24(3), 198–205.

<https://doi.org/10.1097/01.AUD.0000069226.22983.80>

Quar, T.K., Ching, T.Y.C., Newall, P. & Sharma, M. (2013). Evaluation of real-world preferences and performance of hearing aids fitted according to the NAL-NL1 and DSL v5 procedures in children with moderately severe to profound hearing loss. *International Journal of Audiology*, 52: 322-332.

Roush, P. & Jones, C. (2018). Finding the right fit: Pediatric hearing aid coupling options for children. *Pediatric Focus 2*. Accessed from:

https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/resources/evidence/white_paper/documents/technical_paper/Focus_BtB_Pediatric_Pediatric_HA_coupling_options_for_children_210x280_GB_V1.00_028-1904-02.pdf

Scollie, S. (2018). 20Q: Using the Aided Speech Intelligibility Index in Hearing Aid Fittings. *AudiologyOnline*. Article 23707. Accessed from: <https://www.audiologyonline.com/articles/20q-aided-speech-intelligibility-index-23707>.

Scollie, S., Ching, T. Y. C., Seewald, R., Dillon, H., Britton, L., Steinberg, J., & Corcoran, J. (2010). Evaluation of the NAL-NL1 and DSL v4.1 prescriptions for children: Preference in real world use. *International Journal of Audiology*, 49(SUPPL. 1). <https://doi.org/10.3109/14992020903148038>

- Scollie, S., Glista, D., Seto, J., Dunn, A., Schuett, B., Hawkins, M., ... Parsa, V. (2016a). Fitting frequency-lowering signal processing applying the American Academy of audiology pediatric amplification guideline: Updates and protocols. *Journal of the American Academy of Audiology*, 27(3). <https://doi.org/10.3766/jaaa.15059>
- Scollie, S., Levy, C., Pourmand, N., Abbasalipour, P., Bagatto, M., Richert, F., ... Parsa, V. (2016b). Fitting noise management signal processing applying the American Academy of audiology pediatric amplification guideline: Verification protocols. *Journal of the American Academy of Audiology*, 27(3). <https://doi.org/10.3766/jaaa.15060>
- Scollie, S., Seewald, R., Cornelisse, L., Moodie, S., Bagatto, M., Lurnagaray, D., ... Pumford, J. (2005). The desired sensation level multistage input/output algorithm. *Trends in Amplification*, 9(4). <https://doi.org/10.1177/108471380500900403>
- Seewald, R. C., Moodie, K. S., Sinclair, S. T., & Scollie, S. D. (1999). Predictive Validity of a Procedure for Pediatric Hearing Instrument Fitting. *American Journal of Audiology*, 8(2).
- Stelmachowicz, P. G., Pittman, A. L., Hoover, B. M., Lewis, D. E., & Moeller, M. P. (2004). The Importance of High-Frequency Audibility in the Speech and Language Development of Children With Hearing Loss, 130(May 2004), 556–562.
- Uhler, K., Warner-Czyz, A., Gifford, R., & Working Group, P. (2017). Pediatric Minimum Speech Test Battery. *Journal of the American Academy of Audiology*, 28(3), 232–247. <https://doi.org/10.3766/jaaa.15123>
- Walker, E. A., Holte, L., McCreery, R. W., Spratford, M., Page, T., & Moeller, M. P. (2015a). The Influence of Hearing Aid Use on Outcomes of Children With Mild Hearing Loss. *Journal of Speech, Language, and Hearing Research : JSLHR*, 58(5), 1611–1625. https://doi.org/10.1044/2015_JSLHR-H-15-0043
- Walker, E. a., McCreery, R. W., Spratford, M., Oleson, J. J., Van Buren, J., Bentler, R., ... Moeller, M. P. (2015b). Trends and Predictors of Longitudinal Hearing Aid Use for Children Who Are Hard of Hearing. *Ear and Hearing*, 36, 38S–47S. <https://doi.org/10.1097/AUD.0000000000000208>
- Walker, E. A., Curran, M., Spratford, M., & Roush, P. (2019). Remote microphone systems for preschool-age children who are hard of hearing: Access and utilization. *International Journal of Audiology*, 58(4), 200–207. <https://doi.org/10.1080/14992027.2018.1537523>
- Wiesner, Th., Bohnert, A., Limberger, A., Massinger, C., Nickisch, A. (2019, draft). Konsenspapier der DGPP zur Hörgeräte-Versorgung bei Kindern, Vers. 4.0 (Germany) <http://www.dgpp.de/cms/pages/de/profibereich/konsensus.php>
- Wiesner, Th., Boéchat, E. (Brazil), Bohnert, A. (Germany), Chapchap, M. (Brazil), Enderle, A. (Germany), Delaroche, M. (France), Demanez, J.P. (Belgium) , Demanez, L.(Belgium), Gilain, C. (Belgium), Van der Heyden, C. (Belgium), Juarez Sanchez, A. (Spain), Kerkhofs, K. (Belgium), Kerouedan, A. (France),Klinck, V. (Belgium), Leflere, V. (Belgium), Lhussier, TH. (Belgium), Matha , N.(France), Melis, N. (France), Verheyden, P. (Belgium),. Zajicek, F. (Austria). (2018). BIAP Recommendation 12-8. Audiometric procedures in the first year of life Part:12-8.1.4: Auditory brainstem response. Retrieved from: <http://www.biap.org/fr/recommandations/recommandations/tc-12-newborn-hearing-screening-unhs/396-rec-12-8-1-4-en/file>
- Van Eeckhoutte, M., Scollie, S., O'Hagan, R., Glista, D. (2020). Perceptual Benefits of Extended Bandwidth Hearing Aids With Children: A Within-Subject Design Using Clinically Available Hearing Aids. *Journal of Speech, Language, and Hearing Research*, 1-13. https://doi.org/10.1044/2020_JSLHR-20-00271

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