



Roger for home and social situations

Bridging the understanding gap

Introducing Roger

The theory is simple: the more spoken words you can understand, the better you can enjoy the communication.

Roger enables people with hearing difficulties to hear and understand many more words than ever before.

With Roger, you don't need to struggle in difficult situations like noisy restaurants, outdoors or at parties. You can simply listen and live your life to the fullest.



How does it work?

Roger systems feature discreet wireless microphones, used by the person speaking, and tiny Roger receivers that simply click onto your existing hearing aids, cochlear implants or Baha. The result? You hear the speaker's words directly in your ears, without any distracting background noise.



Full speech clarity

Roger digital technology has been scientifically proven to help hearing aid users understand up to 62% more speech in noise and over distance than people without hearing loss.¹


Hassle-free

Roger devices are simple to use. They automatically adjust their own settings to suit the noise and speakers around you.

Highly discreet

Roger devices are small, lightweight and designed with discretion in mind.

¹ Professor Thibodeau, Linda, PhD (2014), Comparison of speech recognition with adaptive digital and FM wireless technology by listeners who use hearing aids, University of Texas, Dallas, USA, The American Journal of Audiology (in press)



Roger microphones



Choose the device that best fits how you live today.



Roger Pen

The Roger Pen is a versatile and stylish wireless microphone that allows you to understand in noise and over distance. It reduces background noise and keeps you connected in one-on-one or group conversations in virtually every listening situation. It features Bluetooth™ for easy cell phone calling.

Roger EasyPen

The Roger EasyPen is designed for you to interact and understand in challenging listening situations. Packed in a stylish design, it features a fully automatic wireless microphone which is easy to use in one-on-one or group conversations.

Roger Clip-On Mic

The Roger Clip-On Mic is a discreet wireless microphone that your conversation partner wears on a shirt, blouse or jacket. Simply clip it on and let the microphone do the rest to help you understand conversations.

Roger receivers



With three types of receivers that function perfectly with hearing aids or implants, there is a Roger receiver available for everyone.



Roger design-integrated receivers

These are designed to click onto Phonak hearing aids. They are also available for selected cochlear implants from Advanced Bionics and Cochlear.



Roger X

This miniature universal Roger receiver is compatible with virtually every behind-the-ear hearing aid, cochlear implant speech processor and can be connected to most Bluetooth streamers.



Roger MyLink

This Roger receiver works with any hearing aid or cochlear implant that features a telecoil.



Easy to use and more flexibility

The Roger Pen and Roger EasyPen are easy to use. They automatically analyze the surrounding noise level and detect their physical position, before using this data to configure their own settings and provide the best speech clarity possible — and all without you pressing a button.

Conference style

Placing on a flat surface, such as a table, the Roger Pen or Roger EasyPen automatically activate the 360-degree (omnidirectional) microphone setting, which picks up voices from all around.

Interview style

When you're in a large group of people, for example at a party, hold the Roger Pen or Roger EasyPen like a reporter's microphone and point it at the person speaking to "zoom in" on their voice.

Around the neck

When it's noisy and you want to focus on the words of just one friend, have your friend hang the Roger Pen or Roger EasyPen around their neck. This focuses the microphone automatically on only their voice (this is also how the Roger Clip-On Mic works).

Situations and recommendations

Learn which Roger microphone best suits each different listening situation.

Situation	Roger Pen	Roger EasyPen	Roger Clip-On Mic
One-on-one conversations – Riding in the car with a friend or shopping	•	•	•
Group conversations* – Noisy restaurants, parties or social gatherings	•	•	
Cell phone calls	•		
Watching TV	•	•	•
Listening to music/multimedia	•	•	•

*Use several microphones at the same time.

One-on-one conversations

Whether at home, in the car, or in public places, Roger technology can help you enjoy crystal clear conversations, without the stress and confusion that background noise brings.

At home

Roger reduces distracting background noise such as the noise of kids playing or the blare of the TV, allowing you to enjoy conversations with family members.

In the car

With engine sounds, road noise and speakers facing away from you, the car can be a tricky listening environment. By bringing the speaker's words directly into your ears, however, Roger overcomes these communication challenges.

In public places

At the mall, walking around town...these situations are full of competing noise which can make it difficult to understand your partner. With Roger, you can capture many more of their words.

Recommended:

For maximum speech clarity, hang the Roger Pen or Roger EasyPen around your friend's neck or clip a Roger Clip-On Mic onto their shirt.



Social situations

Many people with hearing loss struggle to hear and understand in busy social settings. Some even end up avoiding these situations. Roger technology puts the 'life' back into your social life – helping you to connect whatever the background noise.

Dinner with friends or family

Dinner tables and restaurants can often pose a real listening challenge thanks to their complex mix of chatter, background music, and clinking cutlery and glasses. Roger systems reduce these distracting 'other' noises, helping you better catch every word.

At parties

When the music is playing, and people are talking and laughing, it is easier to focus on a conversation with Roger bringing words directly into your ears.

Recommended:

To capture the words of different speakers, either point your Roger Pen or Roger EasyPen at each speaker or place it in the center of the table.



Using technology

Cell phone calls

Enjoy calls again by connecting the Roger Pen to your cell or smartphone via Bluetooth. Or if calling over the Internet, just connect you Roger microphone to your device's headphone output.

Recommended:

Pair the Roger Pen with your Bluetooth cell phone to better catch every word.

Watching TV

Roger microphones can be easily connected to your TV's audio output. This way you can hear your favorite programs without having to turn up the TV.

Recommended:

Use the cable and docking station provided with your Roger microphone to connect it to your TV's headphone jack.

Listening to music

With Roger you can transmit your favorite tunes directly into your hearing aids or CI speech processors. Whatever your preferred music source — smartphone, tablet, PC, Mac® or MP3 player — simply connect your Roger microphone to its headphone output and hear your favorite music.

Using GPS

Never miss a turn again. Connect any Roger microphone to your car's GPS system or smartphone to hear directions clearly in your ears.



Group discussions in noise

When it's very noisy or there are several speakers moving around, catching every word can be a particularly demanding challenge. Roger technology overcomes this challenge by allowing you to use several microphones at the same time.



Simply connect these microphones to your existing hearing technology as normal, then hand a Roger microphone to each speaker to hear every word.



Join us online

There is a whole host of useful Roger-related web pages and tools online. Why not try them today?

HearingLikeMe

An online community for parents, children, teens and anyone whose life has been touched by hearing loss.

www.hearinglikeme.com

Facebook

The very latest Phonak updates in one handy, shareable place.

www.facebook.com/Phonak

Twitter

Bite-sized Phonak updates in 144 characters or less.

www.twitter.com/phonak

YouTube

The official Phonak YouTube channel is full of inspiring and educational content – from product movies and celebrity ambassadors to those heart-warming moments when an infant hears for the first time.

www.youtube.com/phonakofficial



Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

**Interact freely. Communicate with confidence.
Live without limit. Life is on.**

www.phonak-us.com

