


Time to  
get wireless



CROSLink™

# At the heart of innovation

Phonak specializes in the design, development, production and worldwide distribution of technologically advanced wireless communication and hearing systems.

The combination of expertise in hearing technology and a strong co-operation with the hearing healthcare professionals allows Phonak to make a substantial improvement to the quality of life of individuals with hearing impairment and those close to them.

Today, Phonak offers a complete range of digital hearing instruments, along with high-tech specialty products and complementary wireless communication systems. Employing a collective workforce of over 2 600 people, the Phonak Group is active in over 70 countries around the world, and is one of the three leaders in the hearing instrument industry.

**PHONAK**  
hearing systems

[www.phonak.com](http://www.phonak.com)



Phonak wireless  
communication systems



CROSLink™

**PHONAK**  
hearing systems

CROSLink is brought to you by Phonak, the world leader in wireless communication.

This latest innovation is the key to hearing from all directions when one ear is unaidable. CROSLink is a **universal** link, thanks to compatibility with a huge range of BTE instruments. The **wireless** device brings user comfort and freedom of movement, as well as the ultimate in **cosmetic appeal**. This **easy-to-fit** solution gives wearers easy access to a host of previously inaccessible hearing situations.



# Satisfying a huge demand

all-around sound

CROSLink, the unique new universal solution from Phonak. A considerable number of your monaurally-fitted clients are candidates for a BiCROS fitting. Now you can offer an effective, wireless and universal solution to the numerous individuals with an unaidable ear.

Your clients obtain the very best sound quality from their own hearing instruments, and enjoy the benefits of CROS technology without having to make any compromises.

sound quality

universal

wireless



# CROSLink™ at a glance

## Universal

Compatible with most BTE instruments

## Wireless

Reliable performance without an annoying cable

## Cosmetically appealing

Compact design of receiver and transmitter

## Versatile

Ideal for BiCROS or CROS fittings

## Easy to fit

Just clip and go!

## Functional

Volume control and detachable receiver

# individual

## Made to perform

### Versatile – compatible technology

CROSLink meets BiCROS and CROS needs perfectly by setting the switch on the receiver to two dots (••) or one dot (•) respectively. CROSLink can also be used with ear level cochlear implant processors.


### Easy to fit – clip on and go

CROSLink offers a highly effective, yet simple solution. There is no need for complex head measurements or programming. One click and the compact receiver is attached or detached to the audio shoe of any compatible BTE hearing instrument.

# effective

# easy to fit

# one click!



### Functional – tailor made practicality

The volume control on the transmitter allows personalization to different hearing needs and situations. In addition, CROSLink is detachable, so that wearers can remove it easily when the system is not needed.

# enhanced safety



Welcome to the world  
of universal, wireless CROS





# Unwired for sound

CROSLink enables significantly improved communication when one ear is unaidable. It satisfies CROS needs when the other ear is normal, or BiCROS needs when the better ear is impaired.

CROSLink provides more communication enjoyment and greater participation in a variety of environments, such as:

- in a car with friends
- business meetings
- family gatherings
- meals in restaurants

CROSLink also ensures better sound awareness and enhanced safety, e.g. in traffic situations.

**Designed for superb hearing from all directions!**

new design

confidence



# Designed to please

## **Universal – huge choice**

CROSLink opens the door to individual choice and flexibility. You are now free to fit the ideal, state-of-the-art hearing instrument for your client's best ear while providing the benefits of CROS wireless technology. Alternatively, clients can still enjoy the advantages of CROSLink with their existing devices.

## **Wireless – total freedom**

CROSLink is a practical and cosmetically appealing solution when one ear is unaidable. Wearers gain optimal freedom of movement, without the drawbacks of an unsightly and fragile wire.

## **Cosmetically appealing – lifestyle oriented**

The small, sleek CROSLink transmitter has the attractive design of the Phonak Perseo housing. The MicroLink style click-on receiver is very discreet, comfortable and convenient.

# enjoyment



# freedom of movement

