


# Choose the right Dynamic FM receiver



Hearing instrument	ML12i	ML11i	ML10i	ML9i	MLxi*	MyLink+*	Audio Shoe
Exélia Art				■	■	■	AS9
Micro Exélia Art	■				■	■	AS12
Micro Exélia	■				■	■	AS12
Exélia				■	■	■	AS9
Nios micro	■				■	■	AS12
Versáta				■	■	■	AS9
Micro Versáta	■				■	■	AS12
Certéna				■	■	■	AS9
Micro Certéna	■				■	■	AS12
Naída SuperPower		■			■	■	AS11
Naída UltraPower			■		■	■	AS10
Milo micro	■				■	■	AS12
Milo SuperPower		■			■	■	AS11
Milo UltraPower			■		■	■	AS10
Milo Plus FS dAZ						■	-
Milo Plus ITC/HS dAZ						■	-
Milo Plus micro	■				■	■	AS12
Milo Plus SuperPower		■			■	■	AS11
Milo Plus UltraPower			■		■	■	AS10
OK! Plus M		■				■	AS11
OK! Plus SuperPower		■				■	AS11
OK! Plus UltraPower			■			■	AS10
Audéo V / IX						■	-
Audéo YES							-
Savía Art				■	■	■	AS9
Savía				■	■	■	AS9
microPower V / IX						■	-
Eleva				■	■	■	AS9
Valeo					■	■	AS5
miniValeo						■	-
eXtra				■	■	■	AS9
Una				■	■	■	AS9
MAXX					■	■	AS5
PowerMAXX					■	■	AS5
Solo T+					■	■	AS5

\* MLxi can also be used with the majority of BTE hearing instruments produced by other manufacturers.

\*\*MyLink+ can also be used with BTE and ITE hearing instruments produced by other manufacturers, provided they feature a T-coil.

