

Choose the right FM receiver for Cochlear Implant and Baha Models



Cochlear Implant/Baha Model		MLxi	MicroMLxS	MicroLink Freedom	MicroLink CIS MLxi or MicroMLxS	MyLink**	MLxi Baha
Advanced Bionics	Harmony	■	■			■	
	Auria	■	■			■	
	Clarion CII BTE***				■		
	Platinum BTE***				■		
	Platinum (PSP) Bodyworn***				■		
	Clarion S Serie***				■		
	Clarion 1.2***				■		
Cochlear/Nucleus	Nucleus® Freedom™ (BTE)			■		■	
	Nucleus® Freedom™ (bodyworn controller)	■	■			■	
	ESPrIt 3G BTE *	■	■			■	
	ESPrIt/SPrInt***				■		
Cochlear/Baha	Baha Divino						■
	Baha Intenso						■
	Baha Compact						■
Med-EL	Opus 2	■	■			■	
	DUET 2***				■	■	
	TEMPO+/Opus 1 /DUET 1***				■		
MXM	digisonic Digi SP 'K, bodyworn	■	■				
	digisonic Digi SP BTE	■	■				

* 3G adaptor needed

** The orientation of the T-coil in the speech processor, and interference from electrical appliances nearby, may compromise the excellent sound quality of MyLink.

*** Special cables with color tags are needed to connect the MicroLink CIS to the various CI models – please check on www.phonak.com\FM_configurator