

Choose the right FM receiver



| Hearing instrument | MLxi* | ML12i | ML11i | ML10i | ML9i | MicroMLxS* | ML9S | MyLink** | Audio shoe |
|--------------------|-------|-------|-------|-------|------|------------|------|----------|------------|
| Exélia Art | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| Exélia Art micro | ■ | ■ | | | | ■ | | ■ | AS12 |
| Exélia micro | ■ | ■ | | | | ■ | | ■ | AS12 |
| Exélia | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| Nios micro | ■ | ■ | | | | ■ | | ■ | AS12 |
| Versáta | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| Versáta micro | ■ | ■ | | | | ■ | | ■ | AS12 |
| Certéna | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| Certéna micro | ■ | ■ | | | | ■ | | ■ | AS12 |
| Náida SuperPower | ■ | | ■ | | | ■ | | ■ | AS11 |
| Náida UltraPower | ■ | | | ■ | | ■ | | ■ | AS10 |
| Audéo V / IX | | | | | | | | ■ | - |
| Audéo YES | | | | | | | | | - |
| Savia Art | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| Savia | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| microPower V / IX | | | | | | | | ■ | - |
| Eleva | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| Valeo | ■ | | | | | ■ | | ■ | AS5 |
| miniValeo | | | | | | | | ■ | - |
| eXtra | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| Una | ■ | | | | ■ | ■ | ■ | ■ | AS9 |
| MAXX | ■ | | | | | ■ | | ■ | AS5 |
| PowerMAXX | ■ | | | | | ■ | | ■ | AS5 |
| Solo T+ | ■ | | | | | ■ | | ■ | AS5 |

* MLxi and MicroMLxS can also be used with the majority of BTE hearing instruments, providing their manufacturers offer compatible audio shoes.

** MyLink can also be used with other manufacturers' BTE and ITE hearing instruments, provided these feature a T-coil.