


Cable Guide

Product Information – CORE Products

Cable Size



CableFree



	Size 10a*	Size 312*	Size 13*	BTE (CS 44A)	iCube
Exélia Art					
CIC	■				
ITC		■			
HS		■			■
FS			■		■
micro Petite				■	
BTE incl. micro				■	■
Versáta					
CIC	■				
ITC		■			
HS		■			■
FS			■		■
micro Petite				■	
BTE incl. micro				■	■
Certéna					
CIC	■				
ITC		■			
HS		■			■
FS			■		■
micro Petite				■	
BTE incl. micro				■	■
Milo (Plus)					
CIC	■				
ITC		■			
HS		■			
FS			■		
BTE incl. micro				■	
Naída					
BTE				■	■
Nios micro					
BTE				■	■
CRT					
Audéo MINI				■	
Audéo SMART				■	■
Audéo YES				■	■

V1.00/2009-11/12/na

Cable Guide

Product Information – Other Products

Cable Size



CableFree

	Size 10a*	Size 312*	Size 13*	BTE (CS 44A)	iCube
Savia Art					
CIC	■				
ITC		■			
HS		■			
FS			■		
BTE incl. micro				■	
Eleva					
CIC	■				
ITC		■			
HS		■			
FS			■		
BTE incl. micro				■	
eXtra					
CIC	■				
ITC		■			
HS		■			
FS			■		
BTE incl. micro				■	
microPower					
BTE				■	
Una					
CIC	■				
ITC		■			
HS		■			
FS			■		
BTE incl. micro				■	
CRT					
Audéo				■	