

Hearing delight is all around understanding

Phonak CROS provides a truly satisfying wireless solution for wearers with one-sided hearing. Regardless of whether you have a hearing loss in your better hearing ear or not, Phonak CROS offers a solution.



CROS

When your better ear hears normally (CROS)

The transmitter placed behind or inside your poorer ear, picks up sound and transmits it wirelessly to the normal hearing ear. The hearing in your good ear will remain natural and completely unaffected.



BiCROS

When you have a hearing loss in your better ear (BiCROS)

The transmitter placed behind or inside your poorer ear, picks up sound and transmits it wirelessly to the better hearing ear. At the same time, the receiving hearing aid also serves as an amplifier improving hearing in that ear.

Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

Interact freely. Communicate with confidence. Live without limit. Life is on.


Your hearing-care professional:



www.phonak.com

Phonak CROS

Turning one-sided challenges into all around hearing delight



PHONAK

life is on

Introducing all around hearing delight

The new Phonak CROS is the smallest and most stylish wireless CROS/BiCROS solution designed to meet your individual needs and preferences. Phonak CROS consists of only two parts, a transmitter microphone for the ear that cannot benefit from a hearing aid, and a receiver hearing aid on the better hearing ear. This means that when sound comes from your poorer hearing side it will be wirelessly sent to your better ear. This will enable you to hear sounds regardless of where they come from. Experience the new freedom of all around hearing delight.

Phonak CROS is unique versatility and enhanced esthetics

Phonak CROS can be mixed and matched with any style and performance level. The choice is yours.

Phonak CROS is the smallest wireless CROS system. It is available in an award-winning behind the ear style or as a discreet custom in the ear. Vibrant? Discreet? The Phonak CROS transmitters are available in different colors.

Behind-the-Ear



Mocha Taupe

Mid Cloud Gray

Beach Beige

Chocolate Brown



Midnight Black

Silver Shadow

Ocean Blue

In-the-Ear



Hearing delight without compromise

Your benefits	Features
Simulating natural ear acoustics	
The sound from the CROS transmitter is simulated to have traveled through the whole ear and ear canal	Real Ear Sound
Adapting to changing environments	
Enjoy real time and automatic optimization to changing environments	SoundFlow
Full confidence	
Be alert to an unlikely event of a connection interruption or low battery of your CROS transmitter	Signal alerts
One-touch synchronization	
Synchronizes volume or program settings with one touch either on the CROS transmitter or hearing aid	QuickSync