

## Believe in more hearing

Phonak is known for its long-standing commitment and dedication to providing the very best hearing solutions for children. The flexible new Nios S H2O meets children's growing needs and offers an effortless hearing experience at all stages of their active lives.

### Breakthrough audibility for optimal speech and language development

The unique SoundRecover feature provides consistent and effortless audibility of crucial high pitched speech sounds, such as a mother's voice.

### WaterResistant housing for an active childhood

The water and dust resistant housing of Nios S H2O allows children to actively explore the world around them. The Nios S H2O is perfect for activities, such as playing with friends at the beach, splashing in a puddle, enjoying sports or having fun at bath time.

### Small design and best performance

No matter how small, the size of a hearing instrument should never compromise a child's hearing experience. The new Nios S H2O combines the smallest design with the best possible performance.



### Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

**Interact freely. Communicate with confidence. Live without limit.**  
Life is on.

Ask your hearing-care professional:



[www.phonak.com](http://www.phonak.com)

# Nios S H2O

028-1130-02/IV1.00/2011-09/WU Printed in Switzerland © Phonak AG All rights reserved

Believe in more hearing

**PHONAK** life is on

## Hearing all sounds is essential

It is common for children with hearing loss to have less hearing for high pitched sounds. These sounds carry information essential for understanding speech. Often it is not possible to make these sounds sufficiently loud, even with modern hearing instrument technology. Unique to Phonak, SoundRecover takes the softer, high pitched sounds that your child is unable to hear, and shifts them to a region where they can be easily heard.

### The benefits of SoundRecover

Children may experience the following benefits when wearing hearing instruments with SoundRecover:

- Improved hearing of high pitched sounds and better overall understanding (such as the sounds "s" or "f")
- Better speech production
- Improved hearing of everyday environmental sounds (e.g. bird song, alarms, doorbells)

### No whistling, more sound

Hearing aids that whistle (feedback) are not only annoying and embarrassing, they also interfere with hearing. WhistleBlock effectively eliminates whistling whenever it occurs, so your child can benefit from increased amplification everywhere.

## Access to the world of wireless communication

All children should have the freedom to connect to the modern world of communication. Nios S H2O, with its wireless and hands-free solutions, lets your child explore new horizons.

### Dynamic FM – for difficult listening environments

Dynamic FM provides exceptional clarity when your child needs to communicate in difficult listening environments like lively classrooms or busy public places. It also helps where your child doesn't face the speaker such as in the car or in the stroller.

### ComPilot – connecting to a wireless world

As they grow up, children increasingly use telephones, computers, MP3 players and game consoles. ComPilot is the ideal Bluetooth communication interface that lets your child enjoy these devices wirelessly.

### Phonak TVLink S – easy watching with StereoSound

Children love TV, and this dedicated TV solution, combined with ComPilot, offers them an exceptional listening experience. Phonak TVLink S easily connects your child's hearing aids to the TV via Bluetooth, with StereoSound quality.

## Anticipating a child's future needs

Children's needs and the diversity of their listening environments change as they grow up. Phonak has designed its features with a child's future in mind, allowing flexible adjustments of hearing instruments while ensuring consistent high sound quality.



### UltraZoom

Better hearing in background noise

Directional microphones improve hearing ability in background noise. This technology is most beneficial when your child starts to socialize in more noisy environments.



### SoundFlow

Adapting to changes in the listening environment

As children grow up, their listening environment becomes more diverse. SoundFlow accurately identifies subtle changes in the sound environment of your child and automatically adjusts according to the specific situation your child is in.

For more information on the many features and benefits a Phonak can bring to your child visit [www.phonak.com](http://www.phonak.com)