





HE11  
2 cm<sup>3</sup> coupler data

ANSI / ASA S3.22-2014  
IEC 60118-0 : 2015

Output sound pressure level




Acoustic gain



|                              |                    |        |         |         |
|------------------------------|--------------------|--------|---------|---------|
| Frequency range              | <100 Hz - >5000 Hz |        |         |         |
| Total harmonic distortion    | 500 Hz             | 800 Hz | 1600 Hz | 3200 Hz |
|                              | 2.0%               | 1.0%   | 1.0%    | 1.0%    |
| Battery current              | 2.1 mA             |        |         |         |
| Equivalent input noise level | 19 dBSPL           |        |         |         |


Induction coil sensitivity




HE11 680  
2 cm<sup>3</sup> coupler data

ANSI / ASA S3.22-2014  
IEC 60118-0 : 2015

Output sound pressure level




Acoustic gain



|                              |                    |        |         |         |
|------------------------------|--------------------|--------|---------|---------|
| Frequency range              | <100 Hz - >5000 Hz |        |         |         |
| Total harmonic distortion    | 500 Hz             | 800 Hz | 1600 Hz | 3200 Hz |
|                              | 2.0%               | 1.0%   | 1.0%    | 1.0%    |
| Battery current              | 2.0 mA             |        |         |         |
| Equivalent input noise level | 19 dBSPL           |        |         |         |

Induction coil sensitivity



General test information

- Supply voltage 1.3 V / impedance 6.2 Ω.
- Specific measurement settings are used. RTS adjustment with volume control.
- The device is operating in linear mode.
- Low-level expansion is active.
- All data obtained are measured with Phonak Target measurement settings.

Warnings

- ⚠ This hearing aid has an output sound pressure level that can exceed 132 dB SPL. Special care should be taken when fitting this instrument as there is a risk of impairing the residual hearing of the user.
- ⚠ Changes or modifications to the hearing aid that are not explicitly approved by the manufacturer are not permitted. Such changes may damage the ear or the hearing aid.