

Phonak Virto B-Titanium Order Form

Custom Hearing Instruments



Step 1: Order Details

Account Number: 20

Hearing Clinic:

Delivery Address:

Contact Name: Phone Number:

Email Address:

Date Ordered: Date Required:

Step 2: Patient Information

Last Name:

First Name:

Serial number of previous Phonak custom product:

Step 3: HI Warranty / Return Policy

☒ 3 year warranty Extra 1 year warranty (\$140/device)

Step 4: Performance Level

Virto B90 - Titanium - C260 PHO
Virto B70 - Titanium - C261 PHO
VIP option not available

Step 5: AOV Details

Experience Level: First time user Experienced user Age:
(with wearing instruments)

Audiogram (Required for AOV):

Left: HL

250	500	1000	2000	4000


Hz

Right: HL



--	--	--	--	--

Hz

HL - Hearing Loss



Step 6: Product Options

	Smallest possible	Options in the shaded field <u>will</u> increase device size
Power level:	M - level	P-Level SP-Level
Wax system:	<input checked="" type="checkbox"/> Cerustop (CS)	Wax System Cerustop standard option (CS)
Vent type*:	<input checked="" type="checkbox"/> AOV-Acoustical Optimised Venting (AO)	New AOV-O; Open focus occlusion reduction (OO) Customer specific: left____mm right____mm
Cosmetics: Binaural Order	<input checked="" type="checkbox"/> As small as possible – regardless of symmetry (SX)	New Symmetrical build (SY)
Input option:		Telecoil (TP) - (shell style will be changed to CIC)
Removal line:	<input checked="" type="checkbox"/> (RF)	Extended Removal Line +5mm (RFE)
	<input checked="" type="checkbox"/> EasyPhone with Push Button (DT6) 	MiniControl without Push Button (DT5)  Battery door not interchangeable with this option selected

Step 7: Special Instruction

Special instruction / Request:

Step 8: Select your complimentary cleaning kit

C&C ITE Cleaning Kit D-Dry Electronic Dry Aid Kit

Phonak Virto B-Titanium Order Form

Additional Information

PHONAK
life is on

Tips to complete the order form

Step 1 – Order Details

- Please fill in client details including your order number and all relevant dates and contact information

Step 2: Patient Information

- Please provide the client's name plus serial number of their previous Phonak custom product if applicable.

Step 3: HI Warranty / Return Policy

- Standard 3 year warranty or select an extended 4 year warranty.

Step 4: Performance Level

- Select the desired performance level: B90 Premium, B70 Advanced

Step 5: AOV Details

- Indicate the experience level
- Enter Audiometric data (mandatory for AOV)
- Select left/right side

Step 6: Product Options

- Select the custom product options
- The shaded box options will increase device size
- **AOV** = Acoustically Optimised Venting –
- Phonak will set the correct vent diameter and style in relation to Audiogram / shell geometry / feedback canceller and required low frequency gain – this is the recommended vent type.
- **AOV-O** = Acoustically Optimised Venting - Open
- This specific venting is more open and recommended for **first time users** and clients who are sensitive to "occlusion".

Step 7: Special Instruction

- Please let us know about special wishes

Step 8: Select your complimentary cleaning kit

- Select your preferred cleaning kit

Cosmetic aspect

Human ear canals are not fully symmetrical therefore binaural custom in-the-ear devices cannot be built with the same size and position in the ear.

As an option we can design the two devices to look similar in the ear. As a consequence one device may not be as small as possible. Please request this with the symmetrical build (SY) check box.

Shell Styles

IIC



Without push button


CIC




With Telecoil – device will be larger

Power Level

Virto B-Titanium is designed to allow for the smallest device size possible however selecting a P or SP receiver will increase the device's size. AOV-O to be used for hearing loss below 40dB only (see graph).



Impression taking



Recommendation for taking impressions:

- Use EasyView Otoplast
- Use Dreve Otoform A softX 25 shore A (#2956)
- Always fill the complete concha with impression material



EasyView Otoplast
A Phonak invention

Size S # 054-0769 – 20pcs.

M # 054-0768 – 20pcs.

L # 054-0767 – 20pcs.