

Phonak Custom Earpiece 4.0 Order Form

Items on this order form are only compatible with Marvel platform devices.
For older platforms please use the alternative earmould order form.

Step 1: Order details

Account Number: 20

Hearing Clinic:

Delivery Address:

Contact Name: Phone Number:

Email Address:

Date Ordered: Date Required:

VIP 24 hour service (\$50 ex GST)
VIP not possible for Titanium products

Step 2: Patient information

Name:

Audiogram:

	kHz	0.25	0.5	1	2	4
Left HL						
Right HL						

Step 3: Product, style and options



SlimTip 4.0
L R



cShell 4.0
L R



SlimTip 4.0 Titanium
L R



cShell 4.0 Titanium
L R

Please note if a different configuration is required for each ear, please use two separate forms.

	Connection type				Shell colour				Vent style		Shell attribute		
	S Receiver	M Receiver	P Receiver	UP Receiver*	Clear (21)	Rosa Tint (7)	Beige/Pink (26)	Tan (14)	Cavity	Regular	Hollow	Solid	
SlimTip Acrylic (062-0006-01)				-	☒				☒		☒	-	
SlimTip Titanium (062-0015-01)*				-	Only available in Titanium Gray				☒		☒	-	
SlimTip Silicone (062-0007-01)				-	☒		-	-	-	☒	-	☒	
cShell Acrylic (062-0008-01)							☒		-	☒	-	-	
cShell Titanium (062-0009)*					Only available in Titanium Gray				-	☒	-	-	

*Please note SlimTip Titanium and cShell Titanium products are not included free of charge with hearing aid order. *HF3 filter recommended

SlimTip Acrylic options

Wax system ☒ Cerustop (CS) Ext. receiver tube + wax spring (WE) None

SlimTip Titanium options

Wax system ☒ Cerustop (CS) Extra Retention Surface (ERS)

cShell Titanium options

Cerustop (CS) HF 4 (H4) HF 3 (H3) ERS

SlimTip Silicone options

Wax system ☒ None

cShell Acrylic options (HF3 recommended for UP receivers)

Wax system Cerustop (CS) HF 4 (H4) HF 3 (H3) SmartGuard (SG) Ext. receiver tube + wax spring (WE)
Faceplate colour ☒ Pink (26) Cocoa (22) Brown (28) Transparent (13) Tan (14) Black Transparent (47)

Generic options (all earpieces)

Tube length (0-3) Left: _____ Right: _____

Vent size AOV - let Phonak choose. Audiogram req'd. Occluded 1mm 1.1-1.6mm 1.7-2.2mm 2.3-3.0mm Open

Removal line Transparent (21) Removal line extension + 5mm

Locks Canal Skeleton

Step 4: Preferences

If necessary, may we change the following:

- | | | |
|--|-----|----|
| • Step down power level - if audiogram permits (maintain 6dB Gain Reserve) | Yes | No |
| • cShell will be built with an anterior tube outlet in case of anatomical limitations. (See over for more details) | Yes | No |

Special instructions:

Phonak Custom Earpiece 4.0 Order Form

Additional information

Order details

- Please fill in client / order details
- Enter delivery date required
- Enter Audiogram (required for AOV)
- Select warranty and rush order details

Product, Style and Options

- Select earpiece type
- Select receiver type
- Select colours
- Select tube length (mandatory for cShell)
- Select wax system*
- Select vent
- Select removal line

*All SlimTip receivers will be delivered with a CeruShield wax guard in the receiver spout. The SlimTip can be built with an additional wax protection.

Preferences

- Let us know about special wishes
- Select if we may change receiver level or cShell build style under consideration of the audiogram, if necessary
- Select if BTE or xReceiver should be delivered together with the earpiece

Alternative cShell build style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations.



Fig.1 cShell - standard tube





Fig.2 cShell - anterior tube outlet


Tube length

Measure the tube length to select the correct Receiver. Measure both ears for binaural fitting.




Vent Styles

We offer two vent solutions for SlimTip.



Regular:

- Classic vent solution
- Fixed diameter (not enlargable)



Cavity:

- Optimised for best fitrate
- Easy to enlarge by the fitter

Hearing Instrument Compatibility

The SDS 4.0 earpiece portfolio is compatible to Marvel hearing aids or newer. It is not backward compatible to Belong or earlier generation hearing aids.