

Phonak CROS™ B Rechargeable model

- Offers 24 hours* of hearing with a simple 3-hour charge
- The quickest charging, longest lasting rechargeable hearing aid on the market
- Short charge of 30 minutes gives you up to 6 hours of full performance
- Easy handling, no more hassles of disposable batteries
- Groundbreaking built-in rechargeable lithium-ion battery



CROS B-R

* Expected results when fully charged, and up to 80 minutes wireless streaming time. Please refer to www.phonakpro.ca/evidence

Easy-to-use smart charging options

Phonak Charger Case

- A charger, drying kit and protective hard case all in one
- Includes a cleaning tool



Phonak Power Pack

- Easily attached to the Phonak Charger Case
- Ideal for short trips where no power source is available
- Provides 7 full day charges when fully charged



Phonak Mini Charger

- Most compact charging option



Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonak.ca



Phonak CROS™ B-R

When you can follow conversations from whatever direction they come from, life is on

When you can enjoy 24 hours* of hearing with one simple charge, life is on

Phonak CROS™ B-R is the first ever CROS hearing aid to feature a specifically designed **built-in lithium-ion rechargeable battery** that provides 24 hours* of hearing with one simple charge.

Phonak CROS B-R rechargeable hearing aids combine the latest high-performance technology with **groundbreaking battery technology**.

No more hassles with disposable batteries that are small and difficult to handle. Now you can enjoy the freedom and confidence of knowing that CROS B-R is conveniently charged and ready to take you through long days... and nights, if necessary.

* Expected results when fully charged, and up to 80 minutes wireless streaming time. Please refer to www.phonakpro.com/evidence

sonova
HEAR THE WORLD

028-1899-02/2018-05/ Printed in Canada © Sonova AG All rights reserved

A Sonova brand

PHONAK
life is on



What makes CROS B-R so smart

Phonak CROS B-R is based on Phonak Belong™ technology. When CROS B-R is fitted with an Audéo B-R hearing aid, it enables you to hear sounds and voices clearly and follow conversations from whatever direction they come from. Using Binaural VoiceStream Technology™ from Phonak, CROS B-R wirelessly transfers the sound from the non-hearing ear to the hearing aid on the better hearing ear. As a result CROS B-R delivers much better speech understanding in background noise.¹

What makes Belong hearing aids so special is a unique operating system called AutoSense OS™. It senses your surroundings and adjusts your hearing aids every step of the way. Just turn them on and the rest is automatic.

It doesn't matter where you are or what you're doing – CROS B-R and Belong hearing aids were created to simplify your life.

¹ Ebbing, S. & Omisore, D. (2015). Phonak CROS II - Improved speech understanding thanks to binaural beamforming. Phonak Field Study News, retrieved from <https://www.phonakpro.com/us/en/resources/information-forms/evidence.html>, accessed June 14, 2017.



Simply switch them on

CROS B-R is simple and ready to use. Unlike bone-anchored hearing aids which require surgery and an extended time for healing, once CROS B-R is fitted with a Belong hearing aid you can simply switch them on and enjoy hearing from both sides.

With CROS B-R, you will immediately feel included in a conversation, regardless of the situation.

Phonak wireless accessories to boost your hearing aids

Phonak offers a range of wireless accessories to boost the performance of your hearing aids if needed.

To enable you to understand even more while watching TV, talking on the phone or listening to music, we have developed wireless solutions to suit your every need. They also allow you to fully participate in noisy situations or when there's distance between you and the speaker, or when there are multiple speakers.

Phonak CROS B-R in different situations



Conversations in noisy places

Your CROS B-R can zoom in on the voice in front of you and reduce background noise. **This allows you to hear more of what you want and less of what you don't.**



Having a conversation while driving

If your non-hearing ear is closest to your travelling companion, it may be impossible to have a conversation. **With CROS B-R you can easily understand their voice** because it will transfer the sound to your better hearing ear.



Awareness of your environment

With CROS B-R your awareness of sounds improves because you can hear from both directions. It makes conversation easier because you don't rely on people always being on your hearing side.

Top tips

Whether you're getting a CROS B-R for the very first time or upgrading from previous CROS aids, here are our top tips for getting the most out of your Phonak CROS B-R.

1. Talk to your hearing care professional (HCP)

Your HCP will help you by making any necessary adjustments to ensure that your settings are comfortable. Don't be shy in going back to them if you feel something is too loud or "not quite right". They want you to be happy so go back and see them if you need to.

2. Give it time

Adjusting to anything new takes time, whether it's a new pair of shoes or a new car. Adjusting to CROS B-R is no different. Everything will sound different at the beginning so be patient with yourself.

3. Keep trying different environments

Try out familiar situations and then start trying out new environments. Is there any place you've been avoiding because of your hearing loss? See how CROS B-R makes a difference to your comfort level and enjoyment.

4. Tell your friends and family

Let people know that you're adjusting to new hearing aids and it may take some time to get used to.

5. Give yourself a break

If you've had single-sided hearing for a long time, it may take your brain some time to adjust. That is normal. Give yourself a break when needed but remember that the more you use your new CROS B-R, the better it will be.

6. Keep notes

It's a whole new world of sound. If there are any situations where you feel the sound was not quite right—make a note and go back to your HCP to get CROS B-R set up in a way that best fits you.