

Phonak Belong Colours


BTE and RIC



Faceplate colours






Shell colours



Virto B-Titanium



Phonak Coupling Options



Audéo™ B-R Charger Options

Phonak Charger Case RIC

- A charger, drying kit and protective hard case all in one
- Includes a cleaning tool
- Designed for two Audéo B-R hearing aids

Part# 0075-0011



Phonak Mini Charger RIC

- Most compact charging option
- Designed for two Audéo B-R hearing aids

Part# 0075-0003



Universal Phonak Power Pack

- For charging Audéo B-R on the go
- Can be easily attached to the bottom of the Phonak Charger Case
- Provides power for 7 full charges of two Audéo B-R hearing aids

Part# 075-3000-20



Bolero™ B-PR Charger Options

Phonak Charger Case BTE

- A charger, drying kit and protective hard case all in one
- Includes a cleaning tool
- Designed for two Bolero B-PR hearing aids

Part# 0075-0013



Phonak Mini Charger BTE

- Most compact charging option
- Designed for two Bolero B-PR hearing aids

Part# 0075-0014



Universal Phonak Power Pack

- For charging Phonak Bolero B-PR on the go
- Can be easily attached to the bottom of the Phonak Charger Case
- Provides power for 7 full charges of two Bolero B-PR hearing aids

Part# 075-3000-20



Wireless Communication Portfolio

Belong Platform



PilotOne II

The discreet remote control

An easy-to-use remote control for Belong hearing aids. It adjusts the volume and program settings with 4 simple buttons: Home, Volume up, Volume down and Program.



ComPilot II

The multi-purpose streamer

An around the neck streamer and remote control that provides high quality sound through Belong hearing aids when talking on the phone, watching TV or listening to music for a longer period of time.



Phonak TV Connector

Wireless connection to your TV (Audéo Direct only)

Phonak Audéo B-Direct easily connects to your TV or stereo system through a new compact multimedia hub called Phonak TV Connector. This plug and play solution turns your hearing aids into wireless headphones so you can enjoy your favourite TV programs and movies in excellent stereo sound quality.



ComPilot Air II

The clip-on streamer

A small on the go clip-on streamer and remote control focused on Bluetooth technology. It connects Belong hearing aids to cell phones or Bluetooth enabled music devices.



RemoteMic

The lightweight wireless microphone

A Bluetooth clip-on microphone that transmits the speakers' voice directly to both Belong hearing aids via ComPilot II or ComPilot Air II over a distance of up to 20 meters (66 ft).



TVLink II

The TV and music interface

The TVLink II is used with Phonak ComPilot II or ComPilot Air II to turn any Belong hearing aids into a wireless TV headset with stereo sound quality and a range of up to 30 meters (100ft).



EasyCall II

Cell phone solution

This out-of-the-box solution connects Belong hearing aids to any Bluetooth enabled cell phone allowing you to hear a call in both ears.



Roger EasyPen

Microphone

The Roger EasyPen is an exciting product for clients who desire ultimate simplicity with full Roger performance and do not need Bluetooth functionality or manual microphone control.



Roger Table Mic

Microphone

Roger Table Mic is a wireless microphone specifically designed to pick up multiple voices in small and large meetings.



Phonak DECT Phone II

The cordless home phone solution

A cordless phone that transmits sound simultaneously to both ears when Belong hearing aids are detected and used as a regular phone for users with normal hearing.



Roger Pen

Microphone


The Roger Pen is a discreet, cutting-edge microphone with fully automated settings. It features Bluetooth for easy cell phone use and can be used to hear audio from your TV and other multimedia devices.



Roger Clip-On-Mic

Microphone

The Roger Clip-On Mic is a lightweight mini microphone that features advanced microphone processing and is very easy to use.



A Sonova brand

Phonak Belong Models

Phonak Audéo™ B Direct **NEW!**

Audéo B Direct

Battery Size: 13



Phonak Rechargeable

Audéo B-R

Rechargeable
Lithium-Ion Battery



Bolero B-PR

Rechargeable
Lithium-Ion Battery



Phonak Audéo™ B

Audéo B 10

Battery Size: 10



Audéo B 312

Battery Size: 312



Audéo B 312T

Battery Size: 312



Audéo B 13

Battery Size: 13



Phonak Bolero™ B

Bolero B-M

Battery Size: 312



Bolero B-P

Battery Size: 13



Bolero B-SP

Battery Size: 13



Phonak Virto™ B **NEW!**

Non-wireless

Phonak Virto B-Titanium

(B90/B70)

Battery Size: 10



Virto B-IIC

Battery Size: 10



Virto B-10 NW O

Battery Size: 10



Virto B-312 NW O

Battery Size: 312



Wireless

Virto B-10 O

Battery Size: 10



Virto B-10

Battery Size: 10



Virto B-312

Battery Size: 312



Virto B-13

Battery Size: 13



Premium B90

All models

For feature details and power levels per model, please consult Target Software

Direct and Rechargeable options available



Feature	Benefit
Calm Situation	Provides listening comfort and clarity in quiet, one on one conversation.
Speech in Noise	Adapts and reduces noise from behind in real time and applies and optimizes the signal-to-noise ratio (SNR) in noisy environments.
Speech in Loud Noise	Automatically helps you focus on a single speaker directly infront of you in an extremely noisy crowd, for eg. a party or restaurant.
Speech in Car	Reduces car and road noise to make listening in the car more enjoyable.
Comfort in Noise	For listening comfort in loud situations when no speech is present.
Comfort in Echo	Reduces disturbing echo while letting you hear and understand speech naturally.
Music	Offers a full and rich listening experience so you can enjoy your favourite music again.

Speech in Wind	Allows you to enjoy conversation even on a windy day.
EchoBlock	Identifies signals blurred by echo/reverberations and removes them, providing you with clearer hearing.
Speech in 360°	Once activated, this automatically tracks and focuses on the dominant speaker, when you cannot face them, for e.g. driving with passengers behind or beside you.
SoundRelax	Cushions sudden loud or unpleasant impulse sounds such as clanging dishes.
WindBlock	Suppresses wind noise so you can enjoy spending time outdoors.
DuoPhone	The phone signal is automatically heard in both ears, therefore also reducing background noise.
QuickSync	Control both your instruments (program and/or volume changes) with one simple touch.
UltraZoom Premium 33 channel adaptive directional mic	In noisy situations this feature zooms in on voices that are facing you. Noise from the side and back is reduced, for e.g. home activities or in meetings.
SoundRecover2	Makes high pitched sounds easier to hear, for e.g. women's voices.
WhistleBlock 20 channels	Elimination of annoying and embarrassing whistling from the hearing aids before it becomes audible to you or people around you.
NoiseBlock 20 channels	Automatically reduces annoying background noise like engines or fans without reducing the voices you want to hear.
Tinnitus Balance	Every Belong model includes a "built-in" feature designed to reduce the annoyance of your tinnitus, especially in quiet environments.

Don't forget your Wireless Communication Portfolio options (Your practice's standard remote and accessory policy will apply.)

Advanced B70

All models

For feature details and power levels per model, please consult Target Software

Direct and Rechargeable options available



Feature	Benefit
Calm Situation	Provides listening comfort and clarity in quiet, one on one conversation.
Speech in Noise	Adapts and reduces noise from behind in real time and applies and optimizes the signal-to-noise ratio (SNR) in noisy environments.
Comfort in Noise	For listening comfort in loud situations when no speech is present.
Music	Offers a full and rich listening experience so you can enjoy your favourite music again.

Speech in Loud Noise	Helps you focus on a single speaker directly infront of you in an extremely noisy crowd, for eg. a party or restaurant.
Real Ear Sound	This feature helps you regain the ability to tell which direction sounds are coming from the front or back (when wearing behind-the-ear hearing aids).
Speech in 360°	Once activated, this automatically tracks and focuses on the dominant speaker, when you cannot face them, for e.g. driving with passengers behind or beside you.
SoundRelax	Cushions sudden loud or unpleasant impulse sounds such as clanging dishes.
WindBlock	Suppresses wind noise so you can enjoy spending time outdoors.
DuoPhone	The phone signal is automatically heard in both ears, therefore also reducing background noise.
QuickSync	Control both your instruments (program and/or volume changes) with one simple touch.
UltraZoom Advanced 12 channel adaptive directional mic	Automatically zooms in on voices in noisy situations. Noise from the side and back is reduced, for e.g. home activities or in meetings.
SoundRecover2	Makes high pitched sounds easier to hear, for e.g. women's voices.
WhistleBlock 16 channels	Elimination of annoying and embarrassing whistling from the hearing aids before it becomes audible to you or people around you.
NoiseBlock 16 channels	Automatically reduces annoying background noise like engines or fans without reducing the voices you want to hear.
Tinnitus Balance	Every Belong model includes a "built-in" feature designed to reduce the annoyance of your tinnitus, especially in quiet environments.

Don't forget your Wireless Communication Portfolio options (Your practice's standard remote and accessory policy will apply.)

Standard B50

All models

For feature details and power levels per model, please consult Target Software

Direct and Rechargeable options available



Feature	Benefit
Calm Situation	Provides listening comfort and clarity in quiet, one on one conversation.
Speech in Noise	Adapts and reduces noise from behind in real time and applies and optimizes the signal-to-noise ratio (SNR) in noisy environments.
Comfort in Noise	For listening comfort in loud situations when no speech is present.

DuoPhone	The phone signal is automatically heard in both ears, therefore also reducing background noise.
QuickSync	Control both your instruments (program and/or volume changes) with one simple touch.
UltraZoom Standard Single channel adaptive directional mic	Automatically zooms in on voices in noisy situations. Noise from the side and back is reduced, for e.g. home activities or in meetings.
SoundRecover2	Makes high pitched sounds easier to hear, for e.g. women's voices.
WhistleBlock 12 channels	Elimination of annoying and embarrassing whistling from the hearing aids before it becomes audible to you or people around you.
NoiseBlock 12 channels	Automatically reduces annoying background noise like engines or fans without reducing the voices you want to hear.
Tinnitus Balance	Every Belong model includes a "built-in" feature designed to reduce the annoyance of your tinnitus, especially in quiet environments.

Don't forget your Wireless Communication Portfolio options (Your practice's standard remote and accessory policy will apply.)

Essential B30

All models

For feature details and power levels per model, please consult Target Software

Feature	Benefit
Calm Situation	Provides listening comfort and clarity in quiet, one on one conversation.
Speech in Noise	Adapts and reduces noise from behind in real time and applies and optimizes the signal-to-noise ratio (SNR) in noisy environments.

UltraZoom Essential Single channel adaptive directional mic	Automatically zooms in on voices in noisy situations. Noise from the side and back is reduced, for e.g. home activities or in meetings.
SoundRecover2	Makes high pitched sounds easier to hear, for e.g. women's voices.
WhistleBlock 8 channels	Elimination of annoying and embarrassing whistling from the hearing aids before it becomes audible to you or people around you.
NoiseBlock 8 channels	Automatically reduces annoying background noise like engines or fans without reducing the voices you want to hear.
Tinnitus Balance	Every Belong model includes a "built-in" feature designed to reduce the annoyance of your tinnitus, especially in quiet environments.

Don't forget your Wireless Communication Portfolio options (Your practice's standard remote and accessory policy will apply.)