

Phonak Belong platform— technical overview



Phonak Belong product portfolio

Performance level	RIC	BTE	Custom
Premium	Phonak Audéo B90	Phonak Bolero B90	Phonak Virto B90-Titanium
Advanced	Phonak Audéo B70	Phonak Bolero B70	Phonak Virto B70-Titanium
Standard	Phonak Audéo B50	Phonak Bolero B50	
Essential	Phonak Audéo B30	Phonak Bolero B30	

Performance levels Phonak CROS B

Available microphone modes for CROS B

	Automatic StereoZoom	Manual StereoZoom	Real Ear Sound	Omnidirectional
Premium (B90)	•	•	•	•
Advanced (B70)		•	•	•
Standard (B50)			•	•
Essential (B30)				•

New innovation in impression taking – EasyView Otoblock

EasyView Otoblock attaches to your otoscope specula and enables you to take better ear impressions. By providing visual feedback of the ear canal during insertion, EasyView Otoblock can be gently placed directly in front of the ear drum, allowing for a deep ear impression.

EasyView Otoblock provides visualization for the deepest ear impression possible¹— you can insert it past the second bend of the ear canal to receive more detailed anatomical information.

EasyView Otoblock is available in 3 sizes that attach to your otoscope speculum.



¹ Sonova Internal Report, 2016. Impression taking usability study.

Performance levels

Phonak Audéo B and Bolero B	Premium B90	Advanced B70	Standard B50	Essential B30***
AutoSense OS				
Calm Situation	•	•	•	•
Speech in Noise	•	•	•	•
Comfort in Noise	•	•	•	
Music	•	•		
Speech in Loud Noise	•			
Speech in Car	•			
Comfort in Echo	•			
Additional programs				
Max. additional programs	5	5	5	5
Speech in Wind	•			
Comfort in Echo	•			
Speech in Loud Noise	•	•		
Speech in 360°	•	•		
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Comfort in Noise	•	•	•	•
Music	•	•	•	•
Acoustic phone	•	•	•	•
Custom program	•	•	•	•
Streaming programs				
Max. streaming programs	4	4	4	4
Bluetooth audio + mic	•	•	•	•
Bluetooth phone / DECT + mic	•	•	•	•
RemoteMic / Roger	•	•	•	•
Audio jack	•	•	•	•
Features				
UltraZoom	•	•	•	•
SNR-Boost	•	•	•	
FlexControl	•	•	•	
FlexVolume	•	•	•	
DuoPhone	•	•	•	
SoundRecover 2	•	•	•	•
User Preference Tuning	•	•	•	•
Real Ear Sound	•	•	•	
Finetuning channels	20	16	12	8
WhistleBlock	•	•	•	•
NoiseBlock	•	•	•	•
WindBlock	•	•		
EchoBlock	•			
SoundRelax	•	•		
QuickSync	•	•	•	•
AOV	•	•	•	•
Tinnitus Balance	•	•	•	•
auto Acclimatization	•	•	•	•
CROS B compatible****	•	•	•	•
RogerReady	•	•	•	
Roger and directional setting	•	•	•	

Phonak Virto B-Titanium	Premium B90	Advanced B70
AutoSense OS		
Calm Situation	•	•
Speech in Noise	•	•
Comfort in Noise	•	•
Music	•	•
Speech in Car	•	
Comfort in Echo	•	
Additional programs		
Max. additional programs	5	5
Comfort in Echo	•	
Speech in Noise	•	•
Calm Situation	•	•
Comfort in Noise	•	•
Music	•	•
Acoustic phone	•	•
Custom program	•	•
Features		
FlexControl	•	•
FlexVolume	•	•
SoundRecover 2	•	•
User Preference Tuning	•	•
Finetuning channels	20	16
WhistleBlock	•	•
NoiseBlock	•	•
EchoBlock	•	
SoundRelax	•	•
AOV & AOV-O	•	•
Tinnitus Balance	•	•
auto Acclimatization	•	•

*** Not available for Phonak Audéo B-R and Bolero-B-PR.

**** CROS B is not compatible with Audéo B-R and Bolero-B-PR.