

Roger™

(Options 02/03)

Receivers colors

Color overview



		Roger 10 AS10	Roger 11 AS11	Roger 14	Roger 15 AS15	Roger 16 AS16	Roger 17	Roger 18 AS18	Roger 19 AS19	Roger 20	Roger 21
Sand Beige	P1										
Amber Beige	P2										
Sandalwood	P3										
Chestnut	P4										
Champagne	P5										
Silver Gray	P6										
Graphite Gray	P7										
Velvet Black	P8										
Ruby	P9										
Petrol	Q1										
Electric Green	Q2										
Caribbean Pirate	Q3										
Dragon Orange	Q4										
Vanity Pink	Q5										
Blue Lagoon Transparent	Q6										
Royal Purple Transparent	Q7										
Precious Pink	T3										
Beige	O1										
Brown	O3										
Light Gray	O4										
Black	O6										
Pure Transparent	I3										
Blue Transparent	I37										
Purple Transparent	I38										
Gray	I23										
Light Pink	I33										
Light Blue	I34										
Taupe	I35										
Chestnut Brown	I50										
Palladium Black	I70										
Chocolate	E5										
Cliff	E6										
Mid Gray	E8										
Glacier	E9										
Dark Cliff Gray	F0										
Sterling	F3										
Brown	L0										
White	L8										
Black	L9										
Beige/Sand	M1										
Charcoal	M2										
Mocha/Brown	T1										
Smoke/Grey	T2										
Alpine White	XN/T7										
Princess Pink	XP										
Lava Red	M6										
Blue Ocean	M7										
Majesty Purple	M8										
Anthracite	V1										
Beige	V2										
Black	V3										
Ebony	V4										
Nordic Grey	V5										
White	V6										

Roger receiver compatibility guide – Hearing aids

		Design-integrated for Phonak hearing aids						Universal					
		Roger 19	Roger 18	Roger 16	Roger 15	Roger 11	Roger 10	Roger X	Audio shoe	Roger MyLink			
Hearing aids	Belong	RIC		•									
			Audéo B-13		•				•	AS18	•		
			Audéo B-312T						• ¹		•		
			Audéo B-312						• ¹				
			Audéo B-10						• ¹				
			Audéo B-R						• ¹				
		BTE											
			Bolero B-M					• ¹			•		
			Bolero B-P		•			•	AS18	•			
			Bolero B-SP		•			•	AS18	•			
			Bolero B-PR					• ¹					
		Custom											
			Virto B-10 NW O									• ²	
			Virto B-10 O					• ¹			• ²		
			Virto B-10					• ¹			• ²		
			Virto B-312					• ¹			• ²		
			Virto B-13					• ¹			• ²		
		Pediatric											
		Phonak Sky V-M					• ¹				•		
		Phonak Sky V-P		•			•	AS18	•				
		Phonak Sky V-SP		•			•	AS18	•				
		Phonak Sky V-UP	•				•	AS19	•				
		Phonak Sky V-RIC		•			•	AS18	•				
	Power												
		Phonak Naida V-SP		•			•	AS18	•				
		Phonak Naida V-UP	•				•	AS19	•				
		Phonak Naida V-RIC		•			•	AS18	•				
	BTE												
		Phonak Bolero V-M					• ¹				•		
		Phonak Bolero V-P		•			•	AS18	•				
		Phonak Bolero V-SP		•			•	AS18	•				
	RIC												
		Audéo V-13		•			•	AS18	•				
		Audéo V-312T					• ¹			•			
		Audéo V-312					• ¹						
		Audéo V-10					• ¹						
	Custom												
		Virto V-10 O					• ¹				• ²		
		Virto V-10					• ¹			• ²			
		Virto V-312					• ¹			• ²			
		Virto V-13					• ¹			• ²			
	Pediatric												
		Phonak Sky Q-M13				•	•	AS15	•				
		Phonak Sky Q-SP					•	AS11	•				
		Phonak Sky Q-UP					•	AS10	•				
		Phonak Sky Q-RIC				•	•	AS15	•				
	Power												
		Phonak Naida Q-SP					•	AS11	•				
		Phonak Naida Q-UP					•	AS10	•				
		Phonak Naida Q-RIC				•	•	AS15	•				
	BTE												
		Phonak Bolero Q-M 312					• ¹						
		Phonak Bolero Q-M 13				•	•	AS15	•				
		Phonak Bolero Q-P					•	AS16	•				
		Phonak Bolero Q-SP				•	•	AS16	•				
	Basic												
		Phonak Baseo Q-M 13					•	AS9	•				
		Phonak Baseo Q-P					•	AS9	•				
		Phonak Baseo Q-SP					•	AS9	•				
	RIC												
		Audéo Q 312T					• ¹			•			
		Audéo Q 312					• ¹						
	Custom												
		Virto Q 10									•		
		Virto Q 312					• ¹			•			
		Virto Q 13					• ¹			•			

1 Roger X must be used with a ComPilot or ComPilot II streamer
2 only with Telecoil

Some older Spice H20 models are compatible with Roger 15

Roger receiver compatibility guide – Implants

		Design-integrated for implants				Universal					
		Roger 14	Roger 17	Roger 20	Roger 21	Roger X	Audio shoe	Roger MyLink			
Implants	AB	Naida CI Q		•				• ¹		•	
		Harmony / Auria						iConnect		•	
		Neptune						Neptune Connect		•	
	Cochlear	Nucleus 7			•			Mini Microphone 2+		•	•
		Nucleus 5	•					Euro accessory adapter		•	•
		Nucleus 6	•					Euro accessory adapter		•	•
		Baha 5						• ²			
		Baha 4						•			•
		Baha BP100 / BP110						•			•
		Freedom									•
		Med EI	Sonnet				•		FM battery pack cover		•
	Rondo							Mini battery pack		•	•
	Opus2							FM battery pack cover		•	•
	Oticon Medical	Neuro One								•	•
		Ponto 3 / Ponto 3 Power / Ponto 3 SuperPower						Streamer		•	•
		Ponto Plus / Plus Power						Streamer		•	•
		Ponto Pro / Ponto Pro Power									• ³

1 Roger X must be used with a ComPilot or ComPilot II streamer

2 Via Cochlear Mini Microphone 2+

3 Requires external Oticon medical Telecoil

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonakpro.com