

The University of Manchester


MANCHESTER 1824

INFANT HEARING AID ASSESSMENT USING CORTICAL AUDITORY EVOKED POTENTIALS

Kevin J Munro¹ & Suzanne C Purdy²

¹School of Psychological Sciences, University of Manchester
²Psychology Department, University of Auckland

PHONAK, STUTTGART
 24th APRIL 2010



The University of Manchester

MANCHESTER 1824

SEQUENTIAL STEPS

1. HEARING SCREEN
 - e.g., using fixed-level ABR
2. ASSESSMENT
 - e.g., freq-specific ABR or ASSR
3. SELECT & FIT AMPLIFICATION
 - based on good-quality & comprehensive data
4. AUDIBILITY OF SPEECH
5. DISCRIMINABILITY OF STIMULI] CAEPs??
6. MONITOR & PREDICT OUTCOME]

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester



MANCHESTER 1824


WHY OBLIGATORY CAEPS?

1. Originate at end of auditory chain so assessing more of hearing system (& probably more likely to correlate with performance-see later)
2. Corticals work well on awake babies: large amplitude response & easily detected with relatively little repetition
3. Stimuli (such as speech) can be long duration so hearing aid has time to react to real life sounds

(Dillon, 2005)

PHONAK, STUTTGART
24th APRIL 2010






The University of Manchester



MANCHESTER 1824


METHODOLOGY

1. Infant awake on parents lap or in high chair
2. Distracters e.g., toys with lights, mirrors, DVD
3. Snap-on leads so easy to have break e.g. feeding
4. Natural speech stimuli from loudspeaker (m, g, t)
5. 50-100 responses per average
6. Inter-stimulus interval ca 1125 ms
7. Artifact reject at +/- 100 μ V
8. Recording window -100 to +500 ms
9. baseline correct, 30 Hz low pass filter

PHONAK, STUTTGART
24th APRIL 2010







The University of Manchester

MANCHESTER

AUTOMATED WAVEFORM DETECTION & DIFFERENTIATION

(Munro et al, Under Revision)

- n=24 adult listeners
- with and without earplugs to simulate HL
- /m/, /g/ & /t/ at 55, 65 & 75 dB SPL

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester

MANCHESTER 1824

DETECTION OF CAEP (Automated Hotelling T² analysis)

	/m/	/t/	/g/
	Normal		
75 dB SPL	24 (100%)	23 (96%)	23 (96%)
65 dB SPL	24 (100%)	23 (96%)	24 (100%)
55 dB SPL	23 (96%)	23 (96%)	23 (96%)
	Simulated conductive impairment		
75 dB SPL	24 (100%)	24 (100%)	24 (100%)
65 dB SPL	21 (88%)	21 (88%)	23 (96%)
55 dB SPL	16 (67%)	12 (50%)	18 (75%)

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester

MANCHESTER 1824

DIFFERENTIATION OF CAEP (Automated MANOVA analysis)

	/m/ v /t/	/m/ v /g/	/t/ v /g/
	Normal		
75 dB SPL	18/23 (79%)	14/23 (61%)	3/22 (14%)
65 dB SPL	18/23 (79%)	14/24 (58%)	3/23 (14%)
55 dB SPL	16/23 (69%)	15/23 (65%)	5/23 (22%)
	Simulated conductive impairment		
75 dB SPL	8/24 (33%)	13/24 (54%)	5/24 (21%)
65 dB SPL	7/20 (35%)	14/20 (70%)	11/21 (52%)
55 dB SPL	4/11 (11%)	8/14 (60%)	5/11 (46%)

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester

MANCHESTER
1824

EVIDENCE SUPPORTING LINK BETWEEN CAEPS & PERFORMANCE

1. hearing aids
2. cochlear implants
3. auditory neuropathy
4. listening training
5. 'central deafness'

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester



MANCHESTER
1824

RELATIONSHIP BETWEEN OBLIGATORY CAEPS & FUNCTIONAL MEASURES IN YOUNG INFANTS

(Golding et al, 2007)

- Function measured using PEACH questionnaire (in good hearing aid users). Involves parents reporting frequency of aural/oral behaviour
- Quality of CAEP on simple grading
 - 0= no response
 - 1= clear response to one stimuli
 - 2= clear response to 2 stimuli
 - 3 = clear responses to three stimuli

PHONAK, STUTTGART
24th APRIL 2010



The University of Manchester

MANCHESTER
1824

SPEECH PERCEPTION & CAEP IN CHILDREN WITH AUDITORY NEUROPATHY

(Rance et al, 2002)

- Speech recognition & CAEP (measured around the same time) from 15 children with AN who wore HAs
- Around 50% had CAEP and some degree of speech recognition (but wide range of performance)
- The remaining 50% had very poor speech recognition and no CAEP

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester

MANCHESTER
1824

POSSIBLE CLINICAL PROCEDURE

1. To demonstrate audibility
 - /m/, /t/ & /g/ at 65 dB SPL, unaided & aided
2. To evaluate adequacy of frequency response
 - If no response to /m/, /t/ & /g/ consider increasing low, high and mid-frequency gain, respectively
3. To gauge degree of audibility of speech sounds
 - If no response at 65 dB SPL, retest at 75 dB SPL
 - If response at 65 dB SPL, retest at 55 dB SPL

(Adapted from Dillon, 2005)

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester

MANCHESTER 1824

EXAMPLES

PHONAK, STUTTGART
24th APRIL 2010

The University of Manchester

MANCHESTER 1824


Girl with severe hearing-impairment Attends Grammar School for the Deaf

Unaided

Aided

PHONAK, STUTTGART
24th APRIL

Courtesy of Tracy Draper, April 2010 (id: 7)



The University of Manchester


MANCHESTER 1824

Male age 14 years CHARGE syndrome: severe LD & visual impairment Tangible reinforcement audiometry: ca 80 dB HL ?? unaided, 60 dB HL?? aided

	/m/	/g/	/t/
Unaided	RESPONSE [Used info to reduce low frequency gain]	NO RESPONSE	NO RESPONSE
Aided	RESPONSE	RESPONSE	RESPONSE

Courtesy of Vanessa Salisbury, Feb' 2010

PHONAK, STUTTGART
24th APRIL 2010



- The University of Manchester
- MANCHESTER 1824
- ### TAKE-HOME MESSAGE
1. Infants are provided with HAs before good quality behavioural information is available
 2. It is possible to measure speech-evoked CAEPs (aided & unaided) in infants
 3. There may be a role for speech-evoked CAEPs in HA assessment
- PHONAK, STUTTGART
24th APRIL 2010