

# Pick a RIC

## Phonak Audéo™ B – catering to individual needs

### Ease of use

#### Phonak Audéo B-R

- Rechargeable battery
- Easy to charge and use
- Binaural VoiceStream Technology™
- CROS B-R compatible

#### Client needs

- Someone who values convenience and no battery hassles
- Someone who relies on hearing aids to always have battery power



### Direct connectivity to TV and phone

#### Phonak Audéo B-Direct

- Direct connectivity to TV
- Direct connectivity to Bluetooth® enabled phones \*
- Hands-free calls
- Easy to connect and use

\*with Bluetooth® 4.2 wireless technology and most older Bluetooth phones

#### Client needs

- Someone who wants an easy TV solution with outstanding sound quality
- Someone who wants a hands-free phone solution
- Someone who wants to control their hearing aids with an app



### Discreetness

#### Phonak Audéo B-10/312

- Smallest design
- In all performance levels
- Binaural VoiceStream Technology™
- CROS B compatible

#### Client needs

- Someone who wants a small and discreet hearing aid
- An easy to use solution that delivers quality speech understanding in everyday listening situations



### More functions

#### Phonak Audéo B-312T/13

- T-Coil
- More control, longer battery life (13)
- Roger compatible
- In all performance levels
- Binaural VoiceStream Technology™
- CROS B compatible

#### Client needs

- Someone who spends time in groups where it's important to hear in VERY noisy situations
- Someone who needs a T-Coil solution in places of worship, etc.

