


# Phonak Charger Case RIC, Mini Charger RIC and Power Pack


- The **Phonak Charger Case RIC** is a charger, drying kit and protective hard case all in one. Additionally, it includes a cleaning tool in the robust and splash-proof housing designed for two rechargeable hearing aids.
- For charging on the go, the **Phonak Power Pack** can be easily attached to the bottom of the Phonak Charger Case. It provides power for 7 full charges of two rechargeable hearing aids.
- The **Phonak Mini Charger RIC** is a compact charger for rechargeable two hearing aids.

## Description


### Charger Case RIC



### Phonak Power Pack for the Charger Case RIC






### Phonak Mini Charger RIC



# Phonak Charger Case RIC, Mini Charger RIC and Power Pack

## General Information

	Phonak Charger Case RIC	Phonak Power Pack	Phonak Mini Charger RIC
			
Dimensions (L x W x H)	101.3 x 76.1 x 49 mm (4 x 3 x 1.9 ")	84.8 x 59 x 22 mm (3.3 x 2.3 x 0.9 ")	43 x 22 x 29.1 mm (1.7 x 0.9 x 1.1 ")
Weight	90 g (3.1 oz)	37 g (1.3 oz)	10 g (0.4 oz)
Color	Champaign	Champaign	Graphite Gray
Operating conditions	5° to +40° Celsius (41° to 104° Fahrenheit)		
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)		
Humidity during transportation	Up to 90% (non condensing)		
Humidity during storage	0% to 70% (non condensing)		
Humidity during operation	0% to 85% (non condensing)		
Atmospheric pressure during transportation and storage	200 hPa to 1500 hPa		
Atmospheric pressure during use	800 hPa to 1500 hPa		
Included Parts	Power supply mini-USB cable (for charging only) Cleaning tool	mini-USB Cable (for charging only)	Power supply mini-USB Cable (for charging only)
External Interfaces	mini-USB port	mini-USB port	mini-USB port
Battery	None	Li-Polymer	None
Charging time	See Phonak Audéo B Product Information	3 hours	See Phonak Audéo B Product Information
Operating time	See Phonak Audéo B-R Product Information	7 chargers for 2 hearing aids	See Phonak Audéo B-R Product Information

## Power supply

Input	Type: DC, stabilized 5 V +/- 5% Max. current: 500 mA The battery cannot be over-charged. It may become warm during charging. Important notice: It is recommended to only use the original Phonak Chargers and Power Pack
-------	---