

Phonak

Charger Case Combi

Technical Data

Dimensions (L x W x H)	101.3 x 76.1 x 49 mm (4 x 3 x 1.9 ")
Weight	90 g (3.1 oz)
Color	Champaign
Operating conditions	5° to +40° Celsius (41° to 104°Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation	Up to 90% (non condensing)
Humidity during storage	0% to 70% (non condensing)
Humidity during operation	0% to 85% (non condensing)
Atmospheric pressure during use	800 hPa to 1500 hPa
Included Parts	Power supply micro-USB cable (for charging only) Cleaning tool
External Interfaces	micro-USB port
Battery	None
Charging time	See Phonak Audéo M Product Information

Power supply

Input	Type: DC, stabilized 5 V +/- 5% Max. current: 500 mA The battery cannot be over-charged. It may become warm during charging. Important notice: It is recommended to only use the original Phonak Chargers and Power Pack
-------	---



A Sonova brand

PHONAK
life is on