

# Phonak

## Fast Fact

### SWORD™ 3.0

#### What is it?

The 2.4 GHz radio chip, Sonova Wireless One Radio Digital (SWORD) 3.0, has been further developed to unlock the potential of multiple protocols:

- Bluetooth® Classic and Low Energy protocols
- Binaural VoiceStream Technology™
- AirStream™ technology
- Roger™ technology



#### Why use it?

SWORD 3.0 offers

- Binaural streaming of music from any compatible Android™ or iOS® device.
- Binaural phone streaming and real hands-free calling without having to touch their phone.
- Improved speech intelligibility with Binaural VoiceStream Technology™ allowing clients to focus on a single voice in challenging listening situations. This means more signal, less noise, and better understanding.
- AirStream™ technology to deliver top rated sound quality streaming from any TV and stereo system.

- RogerDirect™ allows Roger technology to stream directly to Marvel hearing aids and boosts hearing aid performance in noise and over distance.

#### How to use it?

The implementation of SWORD 3.0 into compatible Phonak hearing aids enables:

- Automatic streaming of the full audio bandwidth in real-time and bi-directionally
- Answering or rejecting phone calls with a simple push of the hearing aid button
- A plug and play TV solution that offers automatic connection to the TV for top rated sound quality streaming from any TV and stereo system
- Easy installation of Roger technology.

\* with Bluetooth® 4.2 wireless technology and most older Bluetooth phones  
The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.  
Apple, the Apple logo, iPhone, and iOS are trademarks of Apple Inc., registered in the U.S. and other countries.  
Android is a trademarks of Google Inc.

For more information, go to [www.phonakpro.com](http://www.phonakpro.com) or contact your Phonak representative