







Product description

Product description

			
	Sky M-M Sky M-M Trial M90/M70/M50/M30	Sky M-PR Sky M-PR Trial M90/M70/M50/M30	Sky M-SP Sky M-SP Trial M90/M70/M50/M30
Battery	312	built-in rechargeable	13
Multi-function button	•	•	•
Telecoil	•		•
Bluetooth®	•	•	•
Nano coating	•	•	•
IP rating	IP68	IP68	IP68
Dimensions (L x W x D)	mm 31.3 x 14.1 x 7.3 inch 1.23 x 0.56 x 0.29"	mm 36.7 x 16.4 x 8.5 inch 1.44 x 0.65 x 0.33"	mm 40 x 16.9 x 8.5 inch 1.57 x 0.67 x 0.33"
Weight	g/oz 2.18 / 0.077	g/oz 3.71 / 0.131	g/oz 4.4 / 0.155
Charging times	0% to 100% 0% to 80% 30 min. fitting	3 h 1.5 h 10 min.	
Max. Power Output (2 cc coupler)	<b>SlimTube 4.0</b> dB SPL 122 58	<b>HE11 680</b> dB SPL 128 63	<b>SlimTube 4.0</b> dB SPL 127 60
Frequency range	Hz <100 - 6500	Hz <100 - 8000	<b>HE11 680</b> dB SPL 130 66
Working current	mA 1.5	mA 1.6	<b>HE11</b> dB SPL 139 81
			<b>HE11 680</b> dB SPL 133 75
			<100 - 5000
			<100 - 5000
			2.1
			2.0



Housing and hook colors

Housing

 H0 Beige	 P1 Sand Beige	 P3* Sandalwood	 P4* Chestnut	 P5* Champagne	 P6 Silver Gray	 P8 Velvet Black
 Q2 Electric Green	 Q3 Caribbean Pirate	 T3 Precious Pink	 T7 Alpine White	 M6 Lava Red	 M7 Blue Ocean	 M8 Majesty Purple

Ear hook

 R1 Green	 R2 Yellow	 R3 Pink	 Q6 Blue	 13 Transparent	 Q7 Purple	 R4 Orange
--	---	---	---	--	---	---

\* The colors P3, P4 and P5 are not available for the performance level M30.

Performance levels

	Premium M90	Advanced M70	Standard M50	Essential M30		Premium M90	Advanced M70	Standard M50	Essential M30
<b>AutoSense Sky OS 3.0</b>					<b>Features</b>				
Calm situation	•	•	•	•	UltraZoom	•	•	•	•
Speech in noise	•	•	•	•	SNR-Boost	•	•	•	
Comfort in noise	•	•	•		DuoPhone	•	•	•	
Music	•	•			SoundRecover2	•	•	•	•
Speech in loud noise <sup>3</sup>	•				User Preference Tuning	•	•	•	•
Speech in car	•				Real Ear Sound	•	•	•	
Comfort in echo	•				Fine tuning channels	20	16	12	8
Media speech + mic	•	•	•	•	WhistleBlock	•	•	•	•
Media music + mic	•	•	•	•	NoiseBlock	•	•	•	•
<b>Additional programs</b>					WindBlock	•	•		
Max. additional programs	3	3	3	3	EchoBlock	•			
Comfort in echo	•				SoundRelax	•	•		
Speech in loud noise	•	•			QuickSync	•	•	•	•
Speech in 360°	•	•			Acoustically Optimized Vent	•	•	•	•
Speech in noise	•	•	•	•	auto Acclimatization	•	•	•	•
Calm situation	•	•	•	•	Environmental balance	•	•	•	•
Comfort in noise	•	•	•	•	RogerDirect™	•	•	•	•
Music	•	•	•	•	Roger and directional setting	•	•	•	
Acoustic phone	•	•	•	•					
Custom program	•	•	•	•					
<b>Streaming programs</b>									
RogerDirect + mic	•	•	•	•					
PartnerMic + mic	•	•	•	•					
Phone call + mic	•	•	•	•					

Fitting

	Sky M-M Sky M-M Trial M90/M70/M50/M30	Sky M-PR Sky M-PR Trial M90/M70/M50/M30	Sky M-SP Sky M-SP Trial M90/M70/M50/M30
Phonak Target	6.1 or higher	6.1 or higher	6.2 and higher
Noahlink Wireless	•	•	•
Noahlink			•
HI-PRO and HI-PRO 2			•

Acoustic coupling

<b>Sound outlet</b>	Standard hook	Damped (HE 11 680) transparent and 6 additional pediatric colors (see page 4) for Sky M-M and Sky M-PR Undamped (HE 11) only in transparent color for Sky M-SP
	Tamperproof hook	Mini damped (HE 11 680) transparent and 6 additional pediatric colors (see page-4) for Sky M-M and Sky M-PR Mini undamped (HE 11) only in transparent color for Sky M-SP
	SlimTube 4.0	Available in left and right and five different lengths (00, 0, 1, 2, 3)
<b>Earpiece for SlimTube</b>	Dome	Cap dome
		Open dome available in three sizes (S, M, L)
		Vented dome available in three sizes (S, M, L)
		Power dome available in three sizes (S, M, L)
	SlimTip™	Acryl Silicone Titanium
<b>Earpiece for hook</b>	Earmold	
	Universal eartip	

Phonak Sky Marvel Trial

The new Phonak Sky M Trial hearing aids can be quickly set to the required performance level in order to suit your recommendation. Various levels of technology in one Phonak Sky M Trial hearing aid provides you with many options for providing the hearing solution at the first appointment or bridging the gap during repairs.

Phonak Charger BTE RIC

This charger is suitable for all rechargeable Marvel BTE and RIC hearing aids. It features a holding space for large earmolds, a USB-C charging inlet and an indicator light.



Phonak Sky M Bluetooth® versions

<b>Direct calls/audio streaming</b>
Bluetooth® 4.2 wireless technology and most older Bluetooth phones

Wireless communication portfolio

		Sky M-M Sky M-M Trial M90/M70/M50/M30	Sky M-PR Sky M-PR Trial M90/M70/M50/M30	Sky M-SP Sky M-SP Trial M90/M70/M50/M30
<b>Wireless accessory</b>	TV Connector	•	•	•
	Phonak PartnerMic™	•	•	•
	Phonak RemoteControl	•	•	•
<b>Roger™</b>	Roger microphones	•	•	•

For more information please visit:



www.phonakpro.com/skymarvel



www.phonakpro.com/roger