

Phonak Fast Facts



Phonak Virto V-10 and Virto V-10 O

What are they?



The Virto V-10 and the Virto V-10 O offer wireless capabilities in very small and discreet packages which feature both a push button and an optional volume control. The Virto V-10 is the smallest full performance hearing instrument with Binaural VoiceStream Technology™ and compatibility with the Phonak Wireless Communication Portfolio. The Virto V-10 O is a wireless CIC capable of QuickSync and compatible with the Phonak Wireless Communication Portfolio.

Why use them?

The Virto V-10 and the Virto V-10 O offer high performance in the smallest size possible.

How do I use them?

Fit the Virto V-10 and Virto V-10 O like any other wireless Phonak hearing instrument using the Target™ 4.2 fitting software and iCube II. Available features will be read out and displayed for each hearing instrument in the Target programming once connected.

Feature	Phonak Virto V-10 	Phonak Virto V-10 O 
Venture Technology	•	•
Floating Antenna™	•	•
AutoSense OS	•	•
Binaural VoiceStream Technology™	•	
Speech in Loud Noise (with adaptive StereoZoom)	•	
Speech in Wind	•	
Speech in 360°	•	
DuoPhone	•	
Tinnitus Balance noise generator	•	•
Push button	•	•
Compatible with CROS II Custom	•	
QuickSync	•	•
Compatible with Wireless Communication Portfolio	•	•
Wireless fitting	•	•

For more detailed information visit: www.phonakpro.com
or contact your Phonak Sales Representative