

Phonak Vitus™



The key to confident hearing

A Sonova brand

PHONAK
life is on

The key to confident hearing

Confidence in technology

As the world's leading provider in the realm of hearing technology, Phonak is dedicated to creating solutions that are not only innovative but also life-changing.

Offering the largest portfolio of hearing solutions available on the market, we cater to all degrees of hearing loss and create products to suit you. With Phonak Vitus and Vitus+, you have the choice of two quality hearing aid families that are based on technology that really works.

We have combined our years of expertise and innovation into a hearing aid that suits your lifestyle and provides you with hearing confidence wherever you go.

The choice is yours

Our new portfolio delivers solid performance in a multitude of hearing situations. When you choose Vitus or Vitus+, you can rest assured that you are getting a hearing aid of excellent value combined with:

- Proven Phonak technology
- Quality sound and hearing performance
- Durable – water and dust resistant

Phonak Vitus

Phonak Vitus provides you with three state-of-the-art, compact Behind-The-Ear models that are ideal for mild to profound hearing loss. With Phonak propriety technology, Vitus delivers when you need it most.


Phonak Vitus+

Phonak Vitus+ offers a discreet Receiver-In-Canal, three In-The-Ear and three Behind-The-Ear hearing aids with an enhanced feature set. Vitus+ is easy to use with intuitive controls and functions automatically in challenging noisy environments.

Packed in an appealing design, there's a Vitus hearing aid to suit your needs.

- Behind-The-Ear models
- Receiver-In-Canal model
- In-The-Ear models

Your hearing loss and cosmetic preference will determine the most suitable Vitus for you.



Key features

Proven Phonak technology ensures that you get quality hearing performance. Vitus and Vitus+ are available with some key features to give your hearing confidence a real boost.



Robust and reliable

Created with robust composite material, Vitus hearing aids are reliable for everyday listening situations. Through rain, snow, heat or dust they are built to meet high quality standards.



Personalized for your needs

Get your hearing care professional to adjust your hearing aid to best suit your hearing needs.

Feature comparison

Below is a feature comparison between Vitus and Vitus+ so you and your hearing care professional can find the one that is right for you.

Feature	Benefit	Vitus	Vitus+
WhistleBlock	No annoying whistling	•	•
UltraZoom (Directional microphone option)	Speech intelligibility in noisy environments	•	•
Everyday Automatic	Automatic activation of UltraZoom in noisy environments		•
SoundRecover	Making high-pitched sounds audible again		•
QuickSync	One press controls both hearing aids		•
Auto Acclimatization	Hearing aid automatically adapts over time for comfortable adjustment		•
Water resistant IP68 ^{1,2}	Water and dust resistant	•	•
Roger compatible ²	Improve speech understanding over distance and in loud noise	•	•
Fine-tuning channels	Flexible adjustments to your individual hearing	4	6
Color options	Matching your preferences	2	5

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 meter of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

² Except Vitus+ In-The-Ear models

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonak.com