




Phonak pediatric solutions

Product and age recommendation overview



PHONAK
life is on



Explore the world with the right equipment

A child is not a small adult, therefore their hearing solution cannot be a one-size-fits-all proposition. The needs of children and how they interact with their surrounding support system of family, friends, and school will change over time. Just as a baby grows from needing a bottle for nourishment to being able cook their own meals as a teen, the hearing equipment that gives them the best and most convenient access to speech, language and communication will change too.

Infants and toddlers are just starting to explore the world of sound and their equipment needs to be secure and easy for parents to handle. As a child moves on to school, their listening situations increase in complexity, and understanding can define academic success. In the teen years, gaining independence and autonomy over their hearing equipment means there are more choices for designing a solution that they love and will use.

The international Phonak Pediatric Advisory Board brings age-appropriate recommendations to the fitting software with Phonak Target Junior mode, but there are many choices for rounding out a child's hearing solution package. The following recommendations for age-appropriate equipment packages can be used as a guide.

Phonak Sky™ V styles



Phonak Sky V-M

- Mini BTE
- 312 battery
- Indicator light
- Tamperproof ear hook & battery door*
- Binaural VoiceStream Technology™
- Telecoil



Phonak Sky V-P, SP, or UP

- Traditional BTE
- 13 or 675 battery
- Indicator light
- Tamperproof ear hook & battery door*
- DAI for ear-level Roger receivers
- Binaural VoiceStream Technology™
- Telecoil



Phonak Sky V-RIC

- Receiver-in-Canal
- 13 battery
- 3 xReceiver options
- DAI for ear-level Roger receivers
- Binaural VoiceStream Technology™
- Telecoil

*Complies with IEC standard (60601-2-66)




Infants and toddlers

Sky V-M	✓ With tamperproof battery door & ear hook
Sky V-P	✓ With tamperproof battery door & ear hook
Sky V-SP	✓ With tamperproof battery door & ear hook
Sky V-UP	✓ With tamperproof battery door & ear hook
Sky V-RIC	
Roger Touchscreen Mic	✓ For use at home and pre-school/childcare
Roger Clip-On Mic	✓ For use at home by parents and siblings
Roger Pen	✓ Experienced parents or children requiring Bluetooth® connection
Roger design-integrated receivers	✓ With tamperproof solution
Roger X	Not recommended for children 0-36 months old
Roger MyLink	
Roger Focus	
CROS II	
DECT Phone II	✓
ComPilot II	
ComPilot Air II	
TVLink II	
PilotOne II	✓ Good for parents
EasyCall II	✓

✓ Recommended

Not recommended



School-aged

Teens

✓	✓
✓	✓
✓	✓
✓	✓
✓ Ages 8 years and older	✓
✓ For school and home use	✓
✓ For home and afterschool use	✓ For home use
✓ For home and afterschool use	✓
✓	✓
✓	✓
✓	✓
✓ Ages 6 and older	✓
✓ For children with a non-aidable ear	✓ For children with a non-aidable ear
✓	✓
✓ Ages 8 years and older or when it can be managed independently	✓
✓ Ages 8 years and older or when it can be managed independently	✓
✓	✓
✓	✓
✓	✓

Roger™ microphones



Roger Touchscreen Mic

Automatic microphones
Small Group mode
10 hour use time
Swipe screen technology
Indicator lights



Roger Pen

Discreet
Automatic microphones
Bluetooth
Rechargeable
Docking station



Roger Clip-On Mic

Small
Directional
Rechargeable
Docking station

Roger receivers



Roger design-integrated

Discreet
Tamperproof design available
Available for selected Phonak traditional BTEs and RICs



Roger X

Small
Universal
Works with most BTEs and CIs
Works with most streamers



Roger MyLink

Universal
Induction
Works with hearing device
T-Coil
Audio output with headphones



Roger Focus

SlimTube
Open fit
312 battery
Normal hearing in at least one ear
• Unilateral hearing loss
• Auditory processing disorder
• Autism spectrum disorder

Wireless accessories



Phonak DECT Phone

Landline wireless phone
Streams to both hearing aids
Booster mode
Also for normal hearing



ComPilot II

Binaural streaming
Stereo
Streams music and phone calls via Bluetooth®
Remote control for hearing aids
24 hour streaming time
Rechargeable



ComPilot Air II

Binaural streaming
Streams music and phone calls via Bluetooth®
Small and flexible to clip
Wireless transmission to hearing aids
Remote control
4 hour streaming time



PilotOne II

Remote control
Small
Easy volume & program selection
AAA battery



Phonak TVLink II

Broadcast stereo audio via Bluetooth®
Digital & analog audio inputs
Easy installation



Phonak EasyCall II

Universal with all Bluetooth® mobile phones
Binaural listening
For phone calls only
10 hour talk time
Hardcase for selected mobile phones available

Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

Interact freely. Communicate with confidence.
Live without limit. Life is on.

www.phonakpro.com/pediatric

