

Phonak Fast Facts



Phonak PilotOne II

What is it?

The Phonak PilotOne II is the easy-to-use remote control for clients with Phonak Venture and Phonak Belong™ hearing aids. It adjusts the volume and program settings of the hearing aids with 4 simple buttons: Home, Volume up, Volume down and Program.

Why use it?

The PilotOne II remote control is designed for simplicity and ease of use. The large buttons allow clients to have easy control over their hearing aids. One additional benefit of PilotOne II is the additional microphone attenuation with a long press of the Volume down. This helps to improve the listening comfort of clients when in noisy situations.



How is it used?

The PilotOne II has to be grouped with Phonak Venture and Phonak Belong hearing aids using Phonak Target 5.0 (or newer) at the end of the fitting session. The functionality of the Home button of PilotOne II can be customized to directly jump to a program of choice.

For more detailed information visit www.phonakpro.com/pilotone