

Phonak

RemoteControl

Technical Data

Dimensions (L x W x H)	84 x 32 x 11.5 mm (3.3 x 1.3 x 0.5")
Weight	26.4 g (9.3 oz), battery included
Color	Champagne/black
Operating conditions	0° to +45° Celsius (32° to 113° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during storage	0% to 70% (non condensing)
Humidity during operation	0% to 95% (non condensing)
Included parts	Lanyard, CR 2032 battery
Ingress protection	IP52

Power supply

Battery	CR 2032
Operating time	12 month at average use

Radio information

Antenna type	Monopole antenna
Operation frequency	2.4 GHz – 2.48 GHz
Modulation	GFSK
Radiated power	< 2.8mW



Standards applied

Emission standards

EN 60601-1-2:2015
IEC 60601-1-2:2014
EN 55011:2009+A1
CISPR11:2009/AMD1:2010
CISPR22:1997
CISPR32:2012
ISO 7637-2:2011
CISPR25:2016
EN 55025:2017

Immunity standards

EN 60601-1-2:2015
IEC 60601-1-2:2014
EN 61000-4-2:2009
IEC 61000-4-2:2008
EN 61000-4-3:2006+A1+A2
IEC 61000-4-3:2006+A1+A2
EN 61000-4-4:2012
IEC 61000-4-4:2012
EN 61000-4-5:2014
IEC 61000-4-5:2014
EN 61000-4-6:2014
IEC 61000-4-6:2013
EN 61000-4-8:2010
IEC 61000-4-8:2009
EN 61000-4-11:2004
IEC 61000-4-11:2004
IEC 60601-1 (§ 4.10.2):2005
ISO 7637-2:2011

Legal manufacturer



Sonova AG
Laubisrütistrasse 28
CH-8712 Stäfa
Switzerland
