

Roger™ receiver overview

(Options 02 and 03)

roger



Roger 10
AS10



Roger 11
AS11



Roger 14



Roger 15
AS15



Roger 16
AS16



Roger 17



Roger 18
AS18



Roger 19
AS19

Color overview

Sand Beige	P1								
Amber Beige	P2								
Sandalwood	P3								
Chestnut	P4								
Champagne	P5								
Silver Gray	P6								
Graphite Gray	P7								
Velvet Black	P8								
Ruby	P9								
Petrol	Q1								
Electric Green	Q2								
Caribbean Pirate	Q3								
Dragon Orange	Q4								
Vanity Pink	Q5								
Blue Lagoon Transparent	Q6								
Royal Purple Transparent	Q7								
Precious Pink	T3								
Alpine White	XN/T7								
Princess Pink	XP								
Lava Red	M6								
Blue Ocean	M7								
Majesty Purple	M8								
Beige	01								
Brown	03								
Light Gray	04								
Black	06								
Pure Transparent	13								
Blue Transparent	37								
Purple Transparent	38								
Gray	23								
Light Pink	33								
Light Blue	34								
Taupe	35								
Chestnut Brown	50								
Chocolate	E5								
Cliff	E6								
Mid Gray	E8								
Glacier	E9								
Sterling	F3								
Brown	L0								
White	L8								
Black	L9								
Beige	M1								
Charcoal	M2								
Mocha	T1								
Smoke	T2								

Roger receiver compatibility guide

			Design-integrated for Phonak hearing aids						Universal			Design-integrated for implants	
			Roger 19	Roger 18	Roger 16	Roger 15	Roger 11	Roger 10	Roger X	Audio shoe	Roger MyLink	Roger 14	Roger 17
Hearing aids	Venture	Pediatric	Phonak Sky V-M						• ¹		•		
			Phonak Sky V-P	•					•	AS18	•		
			Phonak Sky V-SP	•					•	AS18	•		
			Phonak Sky V-UP	•					•	AS19	•		
			Phonak Sky V-RIC		•				•	AS18	•		
			Phonak Naida V-SP		•				•	AS18	•		
		Power	Phonak Naida V-UP	•					•	AS19	•		
			Phonak Naida V-RIC		•				•	AS18	•		
			Phonak Bolero V-M						• ¹		•		
		BTE	Phonak Bolero V-P		•				•	AS18	•		
			Phonak Bolero V-SP		•				•	AS18	•		
			Audéo V-13		•				•	AS18	•		
		RIC	Audéo V-312T						• ¹		•		
			Audéo V-312						• ¹				
			Audéo V-10						• ¹				
		Custom	Virto V-10 O						• ¹		•		
			Virto V-10						• ¹		•		
			Virto V-312						• ¹		•		
			Virto V-13						• ¹		•		
	Quest	Pediatric	Phonak Sky Q-M13				•		•	AS15	•		
			Phonak Sky Q-SP					•	•	AS11	•		
			Phonak Sky Q-UP						•	AS10	•		
			Phonak Sky Q-RIC				•		•	AS15	•		
		Power	Phonak Naida Q-SP					•	•	AS11	•		
			Phonak Naida Q-UP						•	AS10	•		
			Phonak Naida Q-RIC				•		•	AS15	•		
			Phonak Bolero Q-M 312						• ¹		•		
		BTE	Phonak Bolero Q-M 13				•		•	AS15	•		
			Phonak Bolero Q-P			•			•	AS16	•		
			Phonak Bolero Q-SP			•			•	AS16	•		
		Basic	Phonak Baseo Q-M 13				•		•	AS15	•		
			Phonak Baseo Q-P			•			•	AS16	•		
			Phonak Baseo Q-SP			•			•	AS16	•		
		RIC	Audéo Q 312T						• ¹		•		
			Audéo Q 312						• ¹				
		Custom	Virto Q 10										
			Virto Q 312						• ¹		•		
			Virto Q 13						• ¹		•		
Implants	AB	Naida CI Q							• ¹		•		•
		Harmony / Auria							•	iConnect	•		
		Neptune									• ²		
	Cochlear	Nucleus 5							•	Euro accessory adapter	•	•	
		Nucleus 6							•	Euro accessory adapter	•	•	
		Baha 5									• ³		
		Baha 4							•		•		
		Baha BP100 / BP110							•		•		
		Freedom									•		
	Med El	Sonnet							•	FM battery pack cover	•		
		Rondo							•	Mini battery pack	•		
		Opus2							•	FM battery pack cover	•		
	Oticon Medical	Ponto Plus / Plus Power							•	Streamer	•		
		Ponto / Ponto Pro / Ponto Pro Power									• ⁴		

1 Roger X must be used with a ComPilot or ComPilot II streamer

2 Via audio cable connected to external audio input

3 Via audio cable connected to Cochlear Mini Microphone

4 Requires external Oticon medical Telecoil

* Some older Spice H20 models are compatible with Roger 15