

Choose the right FM receiver



Hearing instrument	MLxi*	ML12i	ML11i	ML10i	ML9i	MicroMLxS*	ML9S	MyLink**	Audio shoe
Exélia Art	■				■	■	■	■	AS9
Exélia Art micro	■	■				■		■	AS12
Exélia micro	■	■				■		■	AS12
Exélia	■				■	■	■	■	AS9
Nios micro	■	■				■		■	AS12
Versáta	■				■	■	■	■	AS9
Versáta micro	■	■				■		■	AS12
Certéna	■				■	■	■	■	AS9
Certéna micro	■	■				■		■	AS12
Náida SuperPower	■		■			■		■	AS11
Náida UltraPower	■			■		■		■	AS10
Audéo V / IX								■	-
Audéo YES									-
Savía Art	■				■	■	■	■	AS9
Savía	■				■	■	■	■	AS9
microPower V / IX								■	-
Eleva	■				■	■	■	■	AS9
Valeo	■					■		■	AS5
miniValeo								■	-
eXtra	■				■	■	■	■	AS9
Una	■				■	■	■	■	AS9
MAXX	■					■		■	AS5
PowerMAXX	■					■		■	AS5
Solo T+	■					■		■	AS5

* MLxi and MicroMLxS can also be used with the majority of BTE hearing instruments, providing their manufacturers offer compatible audio shoes.

** MyLink can also be used with other manufacturers' BTE and ITE hearing instruments, provided these feature a T-coil.