

FM Transmitter features

INSPIRO
DISCOVER DYNAMIC FM

DynaMic

SmartLink+

ZoomLink+

EasyLink+



inspiro is a Dynamic FM transmitter for teachers and suitable for use in both regular and special schools.

DynaMic is the passaround microphone for the MultiTalker Network with **inspiro**.

SmartLink+ is the FM transmitter for those who want nothing less than the very best.

ZoomLink+ is the FM transmitter for those who know what they want.

EasyLink+ is the FM transmitter that is easiest to operate.

Dynamic Speech Extractor*	■	■	■	■	■
Monitoring*	■				
DataLogging FM	■	■			
ChannelFinder*	■				
SoundCheck	■				
MultiTalker Network	■	■			
Automatic Scan	■		■		
Direct Frequency Synchronization	■		■	■	■
SoftLanding technology			■	■	■
Frequency change	■		■	■	■
Fast battery charging	■	■	■	■	■
Audio input (TV, radio, computer etc.)	■		■	■	■
Link-up to analog telephone network	■		■	■	■
KeyPad lock	■		■		
Choice of microphone settings (Omni, Zoom, SuperZoom)			■	■	
Remote control (selected Phonak hearing instruments)			■		
Bluetooth for use with cellphones and MP3 players			■		
iLapel	■				
iBoom	■		■	■	■
MM8			■	■	■

* In combination with a MLxi, ML9i, ML10i, ML11i, ML12i, iSense or MLxi Baha FM receiver.

PHONAK

life is on