

## Bridging the understanding gap

### What is Roger?

Even the most powerful hearing aids have limitations. The Roger portfolio consists of wireless microphones and receivers. They pick up the speaker's voice and transmit it to receivers connected to your hearing aids, so you can hear speech clearly over distance and in noise.

Roger is a digital wireless standard specifically designed to help boost your hearing performance<sup>1</sup>.

### New groundbreaking MultiBeam Technology

Roger Select features MultiBeam Technology. By utilizing multiple microphones in six directions, speech from a 360° radius is calculated and compared. The direction with the best clarity is automatically selected.

<sup>1</sup>Mülder, H. (2013). Roger and hearing instruments. Excellent speech understanding in high noise levels. Phonak Field Study News, retrieved from [www.phonakpro.com/evidence](http://www.phonakpro.com/evidence), accessed November 14th, 2017.

## Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

[www.phonak.com](http://www.phonak.com)



## Roger™ wireless microphone solutions

028-1821-02/2018-01/JS Printed in XX © 2018 Sonova AG All rights reserved

A Sonova brand

**sonova**  
HEAR THE WORLD

**PHONAK**  
life is on

# When you feel at the heart of the conversation, life is on

Many people with hearing loss would love to hear well and understand speech in busy social settings.

Roger technology puts the 'life' back into your social life – helping you stay connected regardless of the background noise or listening environment.

## Roger microphones



### Roger Select™

A versatile microphone which automatically selects the person who is talking and with Bluetooth® for phone calls.



### Roger Pen™

A handy microphone with portable design and Bluetooth for phone calls.



### Roger EasyPen

A microphone with only two buttons for ultimate simplicity, but with full Roger performance.

## Roger receivers



### Roger design-integrated

Click onto Phonak hearing aids and selected implants.



### Roger X

Universal miniature receiver.



### Roger MyLink

Neck loop receiver for hearing aids or implants with a telecoil.

Situations with background noise or where the sound source is at a distance are ones where Roger truly excels:

### Dinner with friends or family



Dinner tables and restaurants can pose a real listening challenge. Roger Select cuts out distracting noises and helps you to capture every word.

### At parties



When music is playing and people are talking and laughing, it is easier to focus on a conversation with Roger Pen as it delivers the conversation directly into your ears.

### At home



At home Roger reduces background noise like kitchen equipment or the blare from a TV. This allows you to fully enjoy conversations with family members and friends.


Additional situations where Roger can help overcome communication challenges:

### In the car



Engine sound and road noise make conversations challenging. By bringing the speaker's words directly into your ears, Roger helps overcome the communication challenges.

### Phone calls



Enjoy talking on the phone by connecting the Roger Pen or Roger Select to your phone or computer via Bluetooth or cable.

### Watching TV



Roger microphones can easily be connected to your TV's audio output, so you can hear your favorite programs.