

Phonak

Fast Fact

Phonak Sound Delivery System 4.0

What is it?

Phonak Sound Delivery System 4.0 is the new acoustic coupling portfolio for Marvel hearing aids. It refers to all the parts that support the delivery of sound from the hearing aid to the eardrum. The correct acoustic coupling ensures appropriate acoustic performance in addition to wearing comfort for your clients.

The Phonak Sound Delivery System 4.0 portfolio offers

- newly designed receiver portfolio offering greater reliability
- a new M Receiver as the default receiver
- newly designed domes offering better comfort and retention
- a new cap dome for narrow ear canals
- new custom ear pieces including new silicone and titanium SlimTips.



Why use it?

Benefits for you:

- Flexibility with a broad acoustic coupling portfolio
- Reduced maintenance appointments

Benefits for your customers

- Optimized sound quality due to the new_M receiver
- Comfortable fit even for very narrow ear channels
- For every ear the right material with SlimTip (silicone, titanium and acrylic)
- Longer wearing time between maintenance
- Better fit rate with the new SlimTip titanium and silicon

How to use it?

The Phonak Sound Delivery System 4.0 portfolio offers the following components:

Receiver portfolio : S, M (new), P and UP

Dome portfolio: cap (new), open, vented and power

Custom earpieces : SlimTip (titanium, silicone or acrylic); cShell (acrylic)

Please note that Phonak Sound Delivery System 4.0 is only compatible with Marvel hearing aids and is NOT compatible with older devices.

For more information on the Phonak Sound Delivery System, go to www.phonakpro.com or contact your Phonak representative.