




Phonak Bolero M-M (M90/M70/M50/M30)


SlimTube HE 2 cm³ coupler data

ANSI / ASA S3.22-2014
IEC 60118-0 : 2015

Output sound pressure level




Acoustic gain



| | | | | |
|------------------------------|--------------------|--------|---------|---------|
| Frequency range | <100 Hz - >6500 Hz | | | |
| Total harmonic distortion | 500 Hz | 800 Hz | 1600 Hz | 3200 Hz |
| | 1.5% | 1.5% | 2.0% | 1.0% |
| Battery current | 1.5 mA | | | |
| Equivalent input noise level | 19 dB SPL | | | |


Induction coil sensitivity




HE11 680 2 cm³ coupler data

ANSI / ASA S3.22-2014
IEC 60118-0 : 2015

Output sound pressure level




Acoustic gain



| | | | | |
|------------------------------|--------------------|--------|---------|---------|
| Frequency range | <100 Hz - >8000 Hz | | | |
| Total harmonic distortion | 500 Hz | 800 Hz | 1600 Hz | 3200 Hz |
| | 5.0% | 3.0% | 2.0% | 1.0% |
| Battery current | 1.6 mA | | | |
| Equivalent input noise level | 19 dB SPL | | | |

Induction coil sensitivity



General test information

- Supply voltage 1.3 V / impedance 6.2 Ω.
- Specific measurement settings are used. RTS adjustment with volume control
- The device is operating in linear mode
- Low-level expansion is active
- All data obtained are measured with Phonak Target measurement settings

Warning: Changes or modifications to the hearing aid that are not explicitly approved by the manufacturer are not permitted. Such changes may damage the ear or the hearing aid