

By submitting this order to Phonak for creation of a custom product for your client, you agree that Phonak will scan and then dispose of the ear impression, that the scan and other information generated as part of the order fulfilment process is the property of Phonak and the client has been advised that their personal information (including name, date of birth, gender, device sizing and service history) will be retained for warranty claims, product recall and statistical purposes. The information will only be released or disclosed to another person or agency when specifically requested and with your client's consent. If you would like the ear impression returned, please let us know.

Phonak New Zealand Level 2, The Chester Building,
28 The Warehouse Way, Northcote, Auckland 0627
Fax: 09 486 1895. Email: nzorders@phonak.com

1. Order Details

¹Custom orders: Minimum turn time (in-house): 12 days

Clinic Details	Invoice To _____	HAFS Order
	Ship To _____	Same as Invoice
Client Details Name _____		
Order Details	Date Required _____	Routine ¹ Required by ¹ _____
	Hrg. Instr. Fitter _____	
	Order Number _____	Order Date _____
	Age _____	Enable Order # _____

2. Audiogram

Audiogram details required for all earmould & SlimTip orders including 250Hz

Air

Frequency (Hz): 250, 500, 1k, 2k, 4k, 8k

Hearing Level (dB HL): 0, 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 110

R L

Bone

Frequency (Hz): 500, 1k, 2k, 4k

Hearing Level (dB HL): 0, 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 110

R L

3. Hearing Instruments

	PREMIUM		ADVANCED		STANDARD		ESSENTIAL	
	QTY	COLOUR	QTY	COLOUR	QTY	COLOUR	QTY	COLOUR
Naida P	Naida P90-PR	_____	Naida P70-PR	_____	Naida P50-PR	_____	Naida P30-PR	_____
	Naida P90-UP	_____	Naida P70-UP	_____	Naida P50-UP	_____	Naida P30-UP	_____

4. Tubing options

SlimTube SDS 4.0¹

R	00	0	1	2	3
L	00	0	1	2	3

¹PR models only

5. Domes¹

Right				Left			
Cap	Open	Vented	Power	Cap	Open	Vented	Power
Small							
Medium							
Large							

¹PR models only

6. SlimTips

SlimTip¹

	R	L
SlimTip 4.0 Acrylic (Hollow)		
SlimTip 4.0 Silicone (Solid)		

¹PR models only

7. Earmould

Existing hearing instrument model (to ensure selection compatibility)

	R			L				R			L								
STYLE		Canal			Canal Lock		MATERIAL		Semi-Skeleton		TUBING		Standard - Dry Wall		OPTIONS		Pull Cord (acrylic only)		
		R L			R L			Silicone Soft - 70 Shore				Thick - Dry Wall				Paediatric			Full Helix
		Skeleton			Half Shell Carved			Silicone Soft - 25 Shore				3mm Libby horn tubing				Solid Handle			Removal post (acrylic only)
							Earmould colour (silicone only)												

8. Venting

	R	L
AOV ^{1,2} Audiogram required!		
No Vent		
Pressure (1.2mm)		
Small (1.5mm)		
Medium (2mm)		
Large (2.5mm)		
Other	R _____	L _____

¹ AOV, Acoustically Optimised Venting
² AOV is not recommended for children

9. Accessories

BTE RIC Charging Case
TV Connector
Roger Pen iN _____
Roger Select iN _____
Roger TableMic II iN _____
Roger X _____
PartnerMic
RemoteControl

10. Paediatric Fittings

Tamperproof battery Door Lock ¹ Mini Tone Hook	Junior Care Kit - Kids C & C Kit - ITE & BTE (Standard Tube/Earmould) C & C Kit - ITE & BTE (Slimtubing/RIC)
--	--

¹Fitted as standard when DOB indicates child is <36 Months. Only available with Naida P-UP

Special instructions