

Roger™ receivers overview

(Options 02/03)

roger

Color overview



Roger 10
AS10



Roger 11
AS11



Roger 14



Roger 15
AS15



Roger 16
AS16



Roger 17



Roger 18
AS18



Roger 19
AS19



Roger 20*



Roger 21*

Sand Beige	P1										
Amber Beige	P2										
Sandalwood	P3										
Chestnut	P4										
Champagne	P5										
Silver Gray	P6										
Graphite Gray	P7										
Velvet Black	P8										
Ruby	P9										
Petrol	Q1										
Electric Green	Q2										
Caribbean Pirate	Q3										
Dragon Orange	Q4										
Vanity Pink	Q5										
Blue Lagoon Transparent	Q6										
Royal Purple Transparent	Q7										
Precious Pink	T3										
Beige	O1										
Brown	O3										
Light Gray	O4										
Black	O6										
Pure Transparent	I3										
Blue Transparent	I37										
Purple Transparent	I38										
Gray	I23										
Light Pink	I33										
Light Blue	I34										
Taupe	I35										
Chestnut Brown	I50										
Palladium Black	I70										
Chocolate	E5										
Cliff	E6										
Mid Gray	E8										
Glacier	E9										
Dark Cliff Gray	F0										
Sterling	F3										
Brown	L0										
White	L8										
Black	L9										
Beige/Sand	M1										
Charcoal	M2										
Mocha/Brown	T1										
Smoke/Grey	T2										
Alpine White	XN/T7										
Princess Pink	XP										
Lava Red	M6										
Blue Ocean	M7										
Majesty Purple	M8										
Anthracite	V1										
Beige	V2										
Black	V3										
Ebony	V4										
Nordic Grey	V5										
White	V6										

* 02 only

A Sonova brand

PHONAK
life is on

Roger receiver compatibility guide – Hearing aids

Hearing aids			Design-integrated for Phonak hearing aids		Universal		
			Roger 19	Roger 18	Roger X	Audio shoe	Roger MyLink
Belong™	Pediatric	Phonak Sky™ B-M			• ¹		•
		Phonak Sky B-P		•	•	AS18	•
		Phonak Sky B-SP		•	•	AS18	•
		Phonak Sky B-UP	•		•	AS19	•
		Phonak Sky B-RIC		•	•	AS18	•
		Phonak Sky B-PR			• ¹		•
	Power	Phonak Naída™ B-SP		•	•	AS18	•
		Phonak Naída B-UP	•		•	AS19	•
		Phonak Naída B-R RIC			• ¹		
	BTE	Phonak Bolero™ B-M			• ¹		•
		Phonak Bolero B-P		•	•	AS18	•
		Phonak Bolero B-SP		•	•	AS18	•
		Phonak Bolero B-PR			• ¹		•
	RIC	Phonak Audéo™ B-13		•	•	AS18	•
		Phonak Audéo B-312T			• ¹		•
		Phonak Audéo B-312			• ¹		
		Phonak Audéo B-10			• ¹		
		Phonak Audéo B-R			• ¹		
		Phonak Audéo B-Direct					
	Custom	Phonak Virto™ B-10 NW O					
		Phonak Virto B-10 O			• ¹		
		Phonak Virto B-10			• ¹		• ²
		Phonak Virto B-312			• ¹		• ²
		Phonak Virto B-13			• ¹		• ²
Venture™	Pediatric	Phonak Sky V-M			• ¹		•
		Phonak Sky V-P		•	•	AS18	•
		Phonak Sky V-SP		•	•	AS18	•
		Phonak Sky V-UP	•		•	AS19	•
		Phonak Sky V-RIC		•	•	AS18	•
		Phonak Naída V-SP		•	•	AS18	•
	Power	Phonak Naída V-UP	•		•	AS19	•
		Phonak Naída V-RIC		•	•	AS18	•
		Phonak Bolero V-M			• ¹		•
	BTE	Phonak Bolero V-P		•	•	AS18	•
		Phonak Bolero V-SP		•	•	AS18	•
		Phonak Audéo V-13		•	•	AS18	•
	RIC	Phonak Audéo V-312T			• ¹		•
		Phonak Audéo V-312			• ¹		
		Phonak Audéo V-10			• ¹		
		Phonak Virto V-10 O			• ¹		• ²
	Custom	Phonak Virto V-10			• ¹		• ²
		Phonak Virto V-312			• ¹		• ²
		Phonak Virto V-13			• ¹		• ²

1 Roger X must be used with a ComPilot or ComPilot II streamer
2 only with Telecoil

Roger receiver compatibility guide – Implants

Hearing aids	Quest		Design-integrated for Phonak hearing aids				Universal		
			Roger 16	Roger 15	Roger 11	Roger 10	Roger X	Audio shoe	Roger MyLink
	Pediatric	Phonak Sky Q-M13		•			•	AS15	•
		Phonak Sky Q-SP			•		•	AS11	•
		Phonak Sky Q-UP				•	•	AS10	•
		Phonak Sky Q-RIC		•			•	AS15	•
	Power	Phonak Naida Q-SP			•		•	AS11	•
		Phonak Naida Q-UP				•	•	AS10	•
		Phonak Naida Q-RIC		•			•	AS15	•
	BTE	Phonak Bolero Q-M 312					• ¹		
	RIC	Phonak Audéo Q 312T					• ¹		•
		Phonak Audéo Q 312					• ¹		
	Custom	Phonak Virto Q 10							• ²
		Phonak Virto Q 312					• ¹		• ²
		Phonak Virto Q 13					• ¹		• ²

1 Roger X must be used with a ComPilot or ComPilot II streamer
2 only with Telecoil

Implants			Design-integrated for implants				Universal		
			Roger 14	Roger 17	Roger 20	Roger 21	Roger X	Roger X interface	Roger MyLink
	AB	Naida CI Q		•			• ¹		•
		Harmony / Auria					•	iConnect	•
		Neptune					•	Neptune Connect	
	Cochlear	Nucleus 7			•		•	Mini Microphone 2+	•
		Nucleus 5	•				•	Euro accessory adapter	•
		Nucleus 6	•				•	Euro accessory adapter	•
		Baha 5					•	Mini Microphone 2+	
		Baha 4					•		•
		Baha BP100 / BP110					•		•
		Freedom							•
	MED-EL	SONNET				•	•	FM battery pack cover	•
		RONDO					•	Mini battery pack	•
		RONDO 2							•
		OPUS 2					•	FM battery pack cover	•
	Oticon Medical	Ponto 3 / Ponto 3 Power / Ponto 3 SuperPower					•	Oticon Medical Streamer	
		Ponto Plus / Plus Power					•	Oticon Medical Streamer	
		Ponto Pro / Ponto Pro Power							• ²

1 Roger X must be used with a ComPilot or ComPilot II streamer
2 Requires external Oticon medical Telecoil

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonakpro-us.com