

Phonak Audéo™ Marvel



It's not just a great hearing aid.
It's a multifunctional marvel.

Marvel hearing aids express our commitment to improving hearing without limitations

For over 70 years, Phonak has continuously striven to create technology to improve people's lives. With every new Phonak product, we've pushed the boundaries of hearing technology to create solutions that deliver excellent sound quality, are easy-to-use and provide a true sense of well-being.

Our latest creation is the result of years of extensive research and development all coming together in one marvelous piece of technology. We've focused our attention on what you want and expect from a first-class hearing aid – a clear, rich sound experience.

With this philosophy, combined with modern technology, we've developed a multifunctional hearing aid that conveys a love at first sound experience.



Introducing Phonak Audéo Marvel:

- Clear, rich sound
- Connects to smartphones (iOS, Android and Bluetooth®-enabled devices), TV and more
- Rechargeable models
- Smart apps

Love at first sound... each and every day!

- Exceptional sound quality from the first fit¹
- Better speech understanding in noise²
- Reduced listening effort in noise³
- Top-rated streamed sound quality⁴

Clear, rich sound

From the moment you first put on your Marvel hearing aids, you can expect a sound quality second to none. We've accomplished this by taking the proven performance of our technology and improving it in every way so that your Marvel functions as you'd expect, in even the most challenging environments.

This new generation of Phonak technology recognizes and automatically adapts to precisely match more listening situations than ever before.⁵ You can enjoy better speech understanding in noise², less listening effort³ and be amazed at just how well your hearing aids perform in your everyday life.

Redefining connectivity

Cell phone connectivity

Unlike any other hearing aid available, Marvel can connect directly to your cell phone, including iOS, Android and other Bluetooth®-enabled phones, to stream any kind of audio to your ears. With this freedom, you can enjoy true hands-free phone calls with built-in microphones that pick up your voice. Calls can be answered or declined with a simple push of the button, even if your phone is across the room.

TV and multimedia experience

Combined with universal Bluetooth® connectivity, Phonak Audéo Marvel also features proprietary AirStream™ technology to deliver top-rated streamed sound quality⁴ from any TV or stereo system. Enjoy your favorite TV programs and movies with a simple plug and play solution, called Phonak TV Connector. Marvel hearing aids can distinguish between streamed speech and music signals, and automatically adjust to give you an optimized sound quality.

What you can expect:

- Connectivity to cell phones
- Connectivity to a variety of everyday electronics
- Hearing aids used like wireless stereo headphones



Industry-leading rechargeability

Phonak pioneered rechargeability in the hearing industry with the introduction of lithium-ion rechargeable batteries. Marvel hearing aids are also available with this technology, providing you with a full day of hearing including streaming. This means no more hassles with disposable batteries!

Like your cell phone, Phonak Audéo Marvel devices are always ready to go as soon as you put them on. They switch on automatically when removed from the charger and turn off when placed back in the charging case.

What you can expect:

- Ease of use – no more hassles with disposable batteries
- Full day of hearing including streaming
- Super quick charging
- Battery engineered to last 6 years






Smart apps

Using the **myPhonak app** you get access to the Hearing Diary and Remote Support services. Your hearing care professional has the ability to customize your personal listening experience in real-time, so that you always get the best out of your Marvel hearing aids.

The best part is that for follow-up visits you don't have to physically visit your provider.

The **myCall-to-Text app** provides real-time transcriptions of incoming phone calls, allowing you to communicate more effectively.

The **Remote app** allows for program and volume adjustments from the convenience of your smartphone.



What you can expect:

- Keep track of your hearing experience
- Real-time remote video support from your hearing care professional
- Confidence and peace of mind that your personal hearing needs are addressed⁶
- Hear and read your phone calls to confidently communicate

Audéo Marvel comes in 5 models and 9 colors to suit your individual needs



All models are water
and dust resistant IP68

Battery powered models



Audéo M-312



Audéo M-312T



Audéo M-13T

Rechargeable models



Audéo M-R



Audéo M-RT

Available Fall 2019

Colors



Beige
H0



Sand Beige
P1



Sandalwood
P3



Chestnut
P4



Champagne
P5



Silver Gray
P6



Graphite Gray
P7



Velvet Black
P8



Alpine White
T7

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonak-us.com