

Phonak Virto™ B-Titanium

Product information




Phonak Virto B-Titanium is an innovative custom product offering for mild to severe hearing losses in two performance levels.

Powered by the technology of the Phonak Belong™ platform, Phonak Virto B-Titanium hearing aids are the ideal solution for those clients that prefer the best hearing performance in the most cosmetically attractive custom made hearing solution.

For more information visit www.phonakpro-us.com


Product description

	
	Phonak Virto B-Titanium B90/B70
Power level	
M	•
P	•
SP	•
Hardware configuration / options / technology	
Battery size	10
Push button (for program or volume change)	optional
EasyPhone (for acoustic phone)	optional
MiniControl* (for program or volume change)	optional
Telecoil	optional
Microphone technology	Omni
IP rating	IP68 ¹


*only available without push button

Technical data


Power level		M	P	SP
Maximum Power Output (dB SPL)	2cc coupler	109	115	119
	Ear simulator	120	125	129
Maximum gain (dB)	2cc coupler	40	50	60
	Ear simulator	50	60	70
Frequency range (Hz)	2cc coupler	<100 - 7000	<100 - 6700	<100 - 7000
Working current (mA)	2cc coupler	0.95	1.00	0.95
Fitting range				



For mild to moderate hearing loss, all audiometric configurations



For mild to moderately-severe hearing loss, all audiometric configurations



For moderate to severe hearing loss, all audiometric configurations

Colors




¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 meter of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Fitting

Software	Phonak Target 5.4 or higher
Interfaces	Programming cable, NOAHLink, HI-PRO, HI-PRO2

Wireless Communication Portfolio

Wireless accessory	Phonak Virto B–Titanium B90/B70
	 Phonak MiniControl Optional

Performance levels

	Premium B90	Advanced B70	Features	Premium B90	Advanced B70
AutoSense OS™	Premium	Advanced	FlexControl	•	•
Calm Situation	•	•	FlexVolume	•	•
Comfort in Noise	•	•	SoundRecover 2	•	•
Speech in Noise	•	•	User Preference Tuning	•	•
Music	•	•	Finetuning channels	20	16
Speech in Car	•		WhistleBlock	•	•
Comfort in Echo	•		NoiseBlock	•	•
Additional programs			EchoBlock	•	
Max. additional programs	5	5	SoundRelax	•	•
Comfort in Echo	•		AOV & AOV-0	•	•
Calm Situation	•	•	Tinnitus Balance	•	•
Comfort in Noise	•	•	auto Acclimatization	•	•
Speech in Noise	•	•			
Music	•	•			
Acoustic phone	•	•			
Custom program	•	•			

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally.

www.phonakpro-us.com