
Phonak Charger Case RIC, Mini Charger RIC and Power Pack


- The **Phonak Charger Case RIC** is a charger, drying kit and protective hard case all in one. Additionally, it includes a cleaning tool in the robust and splash-proof housing designed for two rechargeable hearing aids.
- For charging on the go, the **Phonak Power Pack** can be easily attached to the bottom of the Phonak Charger Case. It provides power for 7 full charges of two rechatgeable hearing aids.
- The Phonak Mini Charger RIC is a compact charger for rechargeable two hearing aids.

Description

Charger Case RIC



Phonak Power Pack for the Charger Case RIC



Phonak Mini Charger RIC





Phonak Charger Case RIC, Mini Charger RIC and Power Pack

General Information

	Phonak Charger Case RIC	Phonak Power Pack	Phonak Mini Charger RIC
Dimensions (L x W x H)	101.3 x 76.1 x 49 mm (4 x 3 x 1.9 ")	84.8 x 59 x 22 mm (3.3 x 2.3 x 0.9 ")	43 x 22 x 29.1 mm (1.7 x 0.9 x 1.1 ")
Weight	90 g (3.1 oz)	37 g (1.3 oz)	10 g (0.4 oz)
Color	Champaign	Champaign	Graphite Gray
Operating conditions	5° to +40° Celsius (41° to 104°Fahrenheit)		
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)		
Humidity during transportation		Up to 90% (non condensing)	
Humidity during storage		0% to 70% (non condensing)	
Humidity during operation	0% to 85% (non condensing)		
Atmospheric pressure during transportation and starage		200 hPA to 1500 hPa	
Atmospheric pressure during use		800 hPA to 1500 hPa	
Included Parts	Power supply mini-USB cable (for charging only) Cleaning tool	mini-USB Cable (for charging only)	Power supply mini-USB Cable (for charging only)
External Interfaces	mini-USB port	mini-USB port	mini-USB port
Battery	None	Li-Polymer	None
Charging time	See Phonak Audéo B Product Information	3 hours	See Phonak Audéo B Product Information
Operating time	See Phonak Audéo B-R Product Information	7 charges for 2 hearing aids	See Phonak Audéo B-R Product Information
Power supply	Product Information		Product Information
Input	Type: DC, stabilized 5 V +/- 5% Max. current: 500 mA		
	The battery cannot be over-charged. It may become warm during charging. Important notice: It is recommended to only use the original Phonak Chargers and Pov		

