




## Phonak Audéo L-R (L90/L70/L50/L30/L20/Trial)


**S Receiver**  
**2 cm<sup>3</sup> coupler data**

 ANSI / ASA S3.22-2014  
 IEC 60118-0 : 2015

## Output sound pressure level



## Acoustic gain




|                              |                    |        |         |         |
|------------------------------|--------------------|--------|---------|---------|
| Frequency range              | <100 Hz - >8000 Hz |        |         |         |
| Total harmonic distortion    | 500 Hz             | 800 Hz | 1600 Hz | 3200 Hz |
|                              | 1.5%               | 2.0%   | 2.0%    | 1.0%    |
| Expected operating time*     | 18                 | h      |         |         |
| Equivalent input noise level | 19                 | dBSPL  |         |         |


**M Receiver**  
**2 cm<sup>3</sup> coupler data**

 ANSI / ASA S3.22-2014  
 IEC 60118-0 : 2015

## Output sound pressure level



## Acoustic gain



|                              |                    |        |         |         |
|------------------------------|--------------------|--------|---------|---------|
| Frequency range              | <100 Hz - >8000 Hz |        |         |         |
| Total harmonic distortion    | 500 Hz             | 800 Hz | 1600 Hz | 3200 Hz |
|                              | 1.5%               | 2.0%   | 2.0%    | 1.0%    |
| Expected operating time*     | 18                 | h      |         |         |
| Equivalent input noise level | 19                 | dBSPL  |         |         |

## General test information

- Specific measurement settings are used. RTS adjustment with volume control
- The device is operating in linear mode
- Low-level expansion is active
- All data obtained are measured with Phonak Target measurement settings

## Warnings

- ⚠ This hearing instrument has an output sound pressure level that can exceed 132 dB SPL. Special care should be taken when fitting this instrument as there is a risk of impairing the residual hearing of the user.
- ⚠ Changes or modifications to the hearing aid that are not explicitly approved by the manufacturer are not permitted. Such changes may damage the ear or the hearing aid.


\* Expected operating time of the rechargeable battery depends on active features, the use of wireless accessories, hearing loss, battery age and sound environment.

## Phonak Audéo L-R (L90/L70/L50/L30/L20/Trial)


## MAV Receiver

2 cm<sup>3</sup> coupler dataANSI / ASA S3.22-2014  
IEC 60118-0 : 2015

## Output sound pressure level



## Acoustic gain




|                              |                    |        |         |         |
|------------------------------|--------------------|--------|---------|---------|
| Frequency range              | <100 Hz - >8000 Hz |        |         |         |
| Total harmonic distortion    | 500 Hz             | 800 Hz | 1600 Hz | 3200 Hz |
|                              | 1.5%               | 1.5%   | 1.0%    | 1.0%    |
| Expected operating time*     | 18                 | h      |         |         |
| Equivalent input noise level | 19                 | dBSPL  |         |         |


## P Receiver

2 cm<sup>3</sup> coupler dataANSI / ASA S3.22-2014  
IEC 60118-0 : 2015

## Output sound pressure level



## Acoustic gain



|                              |                   |        |         |         |
|------------------------------|-------------------|--------|---------|---------|
| Frequency range              | <100 Hz - 6300 Hz |        |         |         |
| Total harmonic distortion    | 500 Hz            | 800 Hz | 1600 Hz | 3200 Hz |
|                              | 1.0%              | 1.5%   | 1.0%    | 1.0%    |
| Expected operating time*     | 18                | h      |         |         |
| Equivalent input noise level | 19                | dBSPL  |         |         |


\* Expected operating time of the rechargeable battery depends on active features, the use of wireless accessories, hearing loss, battery age and sound environment.

## Phonak Audéo L-R (L90/L70/L50/L30/L20/Trial)


## UP Receiver

2 cm<sup>3</sup> coupler dataANSI / ASA S3.22-2014  
IEC 60118-0 : 2015

## Output sound pressure level



## Acoustic gain



|                              |                   |        |         |         |
|------------------------------|-------------------|--------|---------|---------|
| Frequency range              | <100 Hz - 6100 Hz |        |         |         |
| Total harmonic distortion    | 500 Hz            | 800 Hz | 1600 Hz | 3200 Hz |
|                              | 1.5%              | 1.5%   | 1.0%    | 1.0%    |
| Expected operating time*     | 18                | h      |         |         |
| Equivalent input noise level | 19                | dBSPL  |         |         |

\* Expected operating time of the rechargeable battery depends on active features, the use of wireless accessories, hearing loss, battery age and sound environment.