

Phonak

Fast Facts

myPhonak 4.0 app

What is it?

The myPhonak 4.0 app is an easy-to-use smartphone app that allows patients to get more from their compatible Phonak hearing aids. The app is packed with enhanced user controls and can also connect patients to their hearing care professionals (HCPs) for a virtual initial fitting or when additional support and adjustments are needed.

How to use it?

myPhonak 4.0 can be downloaded from both the App Store® and Google Play™ Store.

Dynamic Noise Cancellation

- Adjusted on the Speech Focus slider from off to mild, moderate, strong

Tap Control can be configured for:




- Accept/end call or streaming
- Pause/resume streaming
- Access smartphone's voice assistant

Remote Support

- Allows you to fit, fine-tune and support your patients' listening needs remotely in real-time.

Hearing Diary 2.0

- Assists with acclimatization by engaging patients in real-life tasks based on their goals.



Why use it?

All these new features and functions contribute to personalized hearing healthcare. Personalizing can lead to improved patient engagement and increased trust between HCPs and end users which can ultimately improve hearing and communication outcomes (Maidment et al. 2019, Maidment and Ferguson, 2018, Denford S, et al. 2014).

App Store is a service mark of Apple Inc.
Google Play is a trademark of Google LLC.

A Sonova brand

PHONAK
life is on