


# Phonak Paradise RLC coupling guide

## Receivers


**Standard (S)**  
47/111 (2cc)



Mild to moderately-severe hearing loss


### Preferred receiver

**Medium (M)**  
51/114 (2cc)




Mild to severe hearing loss

**Power (P)**  
59/122 (2cc)



Mild to severe hearing loss

**UltraPower (UP)**  
71/132 (2cc)  
*Compatible with cShell only*



Moderate to profound hearing loss

## Phonak Audéo™ P



R



RT



312



13T

## Coupling options



Cap Dome



Open Dome



Vented Dome



Power Dome



CeruShield™



SlimTip  
Acrylic



SlimTip  
Silicone



cShell



cShell  
Titanium