


Phonak RIC xReceiver and Coupling Guide

Fitting Range

xReceiver

Hearing Instrument

Coupling Options



Mild to moderately-severe hearing loss

Standard (xS)
46/111 (2cc)



Phonak Audéo™ B



10 312 312T 13 R Direct
NEW!

Phonak Naída™ V




RIC

Phonak Sky™ V



RIC



Mild to severe hearing loss

Power (xP)
57/124 (2cc)



Phonak Audéo™ B



10 312 312T 13 R Direct
NEW!

Phonak Naída™ V




RIC

Phonak Sky™ V



RIC



Moderate to profound hearing loss

UltraPower (xUP)
66/130 (2cc)



Phonak Audéo™ B



312 312T 13 R Direct
NEW!

Phonak Naída™ V



RIC

Phonak Sky™ V



RIC



Phonak BTE Coupling Guide

Model		Ear Hook	SlimTube	Coupling Options		
Phonak Bolero™ B						
						
M	P		HE 10 (unfiltered)			
			HE 10 680 (filtered)			
PR	SP					
Phonak Sky V						
						
M	P					
Phonak Naída V						
			HE 11 (unfiltered)			
SP	UP		HE 11 680 (filtered)	Power SlimTube II (00, 0, 1, 2, 3)	Traditional Earmold*	
Phonak Sky V						
						
SP	UP					
Phonak CROS B						
			CROS Hook (13 models only)			
312	13					
Phonak CROS II						
						
312	13					

* Available from your preferred earmold manufacturer