


Your hearing. Your solution.




Patient Name

Date

Provider Name

Next appointment date/time

Your audiogram



Hearing aid styles

Extended Wear



Custom

(In-The-Ear)



RIC

(Receiver-In-Canal)




BTE

(Behind-The-Ear)



RIC and BTE available in a host of hair, skin, fashion and fun colors:



Select your hearing goals

Phonak has the right hearing solution for you – regardless of your type of hearing loss, lifestyle, personal preferences, age or budget.

Premium



Follow discussions in the car



Enjoy conversations at noisy social events



Engage in outdoor sports activities



Enjoy a visit to the museum

Premium

Developed to satisfy those who desire and expect nothing but the best. Thanks to a complete set of trailblazing solutions that uniquely address all hearing goals, these hearing aids support understanding in even the most difficult listening situations, for example at a busy social event.

Advanced



Hear voices from any direction



Enjoy fuller and richer music



Enjoy conversations at restaurants

Advanced

The perfect balance of performance and price. An impressive offering of features that enhances hearing ability in a variety of listening environments, including those with moderate levels of background noise, for example at a restaurant.

Standard



Converse with colleagues



Chat on the phone

Standard

Ideally suited for people who expect value in state-of-the-art technology. A competitive feature set that has the capability to support the hearing demands of daily life with minimal to moderate levels of background noise, for example using the telephone in a public place.

Essential



Enjoy personal discussions

Essential

The Essential products are designed for people looking for an affordable yet effective solution to enhance daily communication abilities, mostly in quiet or minimal background noise, for example one on one conversations.

Phonak Wireless Communication Portfolio

Control

Easy to use remote control for your hearing aids and an App for discreet access.



PilotOne
PilotOne II



RemoteControl App

Phone

Talk on the phone at home and on your cell phone with ease.



DECT
DECT II



EasyCall
EasyCall II

Connection

Wirelessly stream TV, music, and conversations in stereo quality.



ComPilot Air II



ComPilot
ComPilot II



TVLink II

Noise

Understand speech over distance in groups or one-on-one in the most challenging, noisy environments.



Roger Pen
Roger EasyPen



Roger Clip-On Mic



Roger Receiver

Distance

Enhance one-on-one conversations over distance.



RemoteMic

Phonak feature overview

AutoSense OS

Accurately and precisely adapts to provide a seamless listening experience
Move from one environment to the other, no button-pushing needed

	Premium	Advanced	Standard	Essential
 Speech in Loud Noise (automatic) Automatically focuses on a single voice directly in front of you in a noisy crowd Gather with your family at a restaurant	Yes	No	No	No
 Speech in Car Reduces road and car noise to make conversations in the car easier Enjoy catching up with passengers or a peaceful solo ride	Yes	No	No	No
 Comfort in Echo Automatically reduces reverberations to improve understanding in places with echos Follow the tour guide at the museum without difficulty	Yes	No	No	No
 Speech in Loud Noise (manual) Automatically focuses on a single voice directly in front of you in a noisy crowd Gather with your family at a restaurant	Yes	Yes	No	No
 Music Allows for a more natural music listening experience Enjoy the melodies and harmonies at the symphony	Yes	Yes	No	No
 Comfort in Noise Prevents the hearing aids from becoming too loud in noisy environments Appreciate the cheering crowd at a ball game	Yes	Yes	Yes	No
 Speech in Noise Follow voices in a noisy environment Enjoy conversations at a wedding reception	Yes	Yes	Yes	Yes
 Calm Situations Engage in one on one conversations in quiet Catch up with loved ones in your living room	Yes	Yes	Yes	Yes

Phonak features

 Speech in Wind Allows you to enjoy conversations, even on a windy day Hear the conversation on the golf course	Yes	No	No	No
 Speech in 360° Automatically focuses on the loudest speaker, regardless of their location Follow conversations while driving in the car	Yes	Yes	No	No
 SoundRelax For a more relaxing listening experience Hear beyond clattering dishes or slamming doors	Yes	Yes	No	No
 DuoPhone Hear phone conversations in both ears Talk on your home or cell phone	Yes	Yes	Yes	No
 SoundRecover Gives more high-frequency audibility Enjoy birds chirping, understand children's voices more clearly	Yes	Yes	Yes	Yes
 WhistleBlock Eliminates annoying whistling from hearing instruments Keep focused on the conversation without the worry of interruption	Yes	Yes	Yes	Yes

