

# Custom Ear Piece Order Form

Venture



## Customer Information

Ship To Account Number: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Bill To Account Number: \_\_\_\_\_

Third Party Patient Number: \_\_\_\_\_

Date: \_\_\_\_\_

Purchase Order Number: \_\_\_\_\_

Contact Name: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

## Patient Information

Last Name: \_\_\_\_\_

First Name: \_\_\_\_\_

Age: \_\_\_\_\_

Gender: \_\_\_\_\_

## Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

## HI Warranty/Rush Options

2<sup>nd</sup> year    3<sup>rd</sup> year    4<sup>th</sup> year /    24-hour service (\$59)

## Technology Level

**V90 (Premium)**    **V70 (Advanced)**    **V50 (Standard)**    **V30 (Essential)** N/A for Sky V

## Hearing Device Selection - RIC

Qty	Instrument Type	Instrument Style			
		10	312	312T	13
_____	Phonak Audéo V	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	Phonak Naída V-RIC	N/A	N/A	N/A	<input type="checkbox"/>
_____	Phonak Sky V-RIC	N/A	N/A	N/A	<input type="checkbox"/>

Specify color choice: \_\_\_\_\_

## Custom Ear Piece Selection - RIC

Side	Shell Styles	Power Level			
L	R	xS (46/114)	xP (57/127)	xSP+ (63/131)	xUP (66/130)
<input type="checkbox"/>	<input type="checkbox"/> SlimTip Soft	<input type="checkbox"/>	N/A	N/A	N/A
<input type="checkbox"/>	<input type="checkbox"/> SlimTip Hard	<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A
<input type="checkbox"/>	<input type="checkbox"/> cShell Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Receiver Length (0-3): R \_\_\_\_\_ L \_\_\_\_\_

\*Audeo V-312, 312T and 13 only  
\*\*Naída V-RIC and Sky V-RIC only

## Hearing Device Selection - BTE

Qty	Instrument Type	Instrument Style		
		M	P	SP
_____	Phonak Bolero V	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
_____	Phonak Sky V	<input type="checkbox"/>	<input type="checkbox"/>	N/A

Specify color choice: \_\_\_\_\_

## Custom Ear Piece Selection - BTE

Side	Shell Styles
L	R
<input type="checkbox"/>	<input type="checkbox"/> SlimTip Soft
<input type="checkbox"/>	<input type="checkbox"/> SlimTip Hard

One SlimTube is included at no charge (default: size 2)

SlimTube Length (0-3): R \_\_\_\_\_ L \_\_\_\_\_

## Hearing Device Selection - CROS

Qty	Instrument Type	Instrument Style	
		312	13
_____	Phonak CROS II	<input type="checkbox"/>	<input type="checkbox"/>

Specify color choice: \_\_\_\_\_

## Custom Ear Piece Selection - CROS

Side	Shell Styles
L	R
<input type="checkbox"/>	<input type="checkbox"/> CROSTip Hard

Wire/Tube length (0-3): \_\_\_\_\_

## Colors

P1	P2	P3	P4	P5	P6	P7	P8	P9	Q1	O1
Sand Beige	Amber Beige	Sandalwood	Chestnut	Champagne	Silver Gray	Graphite Gray*	Velvet Black	Ruby*	Petrol	Beige
Exclusive to Phonak Sky V:		Q2	Q3	T3	T7	M6	M7	M8		
		Electric Green	Caribbean Pirate	Precious Pink	Alpine White	Lava Red	Blue Ocean	Majesty Purple	*Not available on Phonak Sky V	

## Ear Piece Options

### Shell Color Options

	Pink		Tan		Cocoa		Brown		Transp.		Other*	
Shell Styles:	L	R	L	R	L	R	L	R	L	R	L	R
SlimTip Soft	N/A		N/A		N/A		N/A		Std.		N/A	
SlimTip Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
cShell Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CROSTip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

\*Other color options: Red, Blue

Specify color choice: \_\_\_\_\_

### Faceplate Color Options

	Pink		Tan		Cocoa		Brown	
Shell Styles:	L	R	L	R	L	R	L	R
cShell Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Shell Options

	Canal Lock*		Skeleton Lock*		Helix Lock*		Extra Retention		Extra Large	
Shell Styles:	L	R	L	R	L	R	L	R	L	R
SlimTip Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A	
cShell Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A		<input type="checkbox"/>	<input type="checkbox"/>
CROSTip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A		N/A	

\*The shell and lock will be the same color

### Venting Options

	AOV (Audiogram required)		Pressure Vent		2.0mm SAV		2.5mm SAV		3.0mm SAV		Customized Large Vent		Cavity Vent		No Vent	
Shell Styles:	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R
SlimTip Soft	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A		<input type="checkbox"/>	<input type="checkbox"/>	N/A		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SlimTip Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
cShell Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Wax Options

	Cerustop		Extended Receiver Tube		Spring Guard*		No Wax Prevention*	
Shell Styles:	L	R	L	R	L	R	L	R
SlimTip Hard (w/ Extra Retention)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
cShell Hard	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

\*xUP cShell Hard only compatible with these wax options.

### Removal Line Color Options

	L	R
Transparent (Std.)	<input type="checkbox"/>	<input type="checkbox"/>
Tan/Pink	<input type="checkbox"/>	<input type="checkbox"/>
Brown/Cocoa	<input type="checkbox"/>	<input type="checkbox"/>
No Removal Line	<input type="checkbox"/>	<input type="checkbox"/>

## Accessories

- Conversation Set (ComPilot Air II/RemoteMic)
- Entertainment Set (ComPilot II/TVLink II)
- ComPilot Air II
- ComPilot II
- TVLink II
- RemoteMic
- PilotOne II
- DECT II

## Preferences

- If necessary, may we change the following:  Please Call
- Yes No
- xUP/xSP plus to xP **OR** xP to xS *if audiogram permits*
- Build larger if required

## Measuring Tips

1. Make sure you are using the Phonak measuring tool
2. Make sure Receiver/SlimTube length measuring tool is parallel to the floor (straight on top of the ear)
3. Do not apply pressure to the tool (let it sit on the ear)
4. Use block colors to determine size (if on the line, use the smaller size)
5. Measure to the top of the ear canal opening

## Special Instructions



All of our products, including custom products and spare parts, can be ordered online in the Phonak store.

Sonova USA Inc. is not responsible and assumes no liability for any non-Phonak manufactured device or accessory sent by you to Phonak. Please ensure that you only include Phonak devices and accessories herein. Products, options and accessories are subject to change without notice.